



PRODUCT	PRODUCT ID	MODULE	CATEGORY	SIZE	TOE CONFIGURATION
Agilix	F15	00	01 - 09	22 - 31	RU: Regular, Universal SR: Sandal, Right* SL: Sandal, Left* <i>*Sandal Toe Available Sizes:</i> 23 - 28: ShockWave 24 - 28: Pacifica & Kinnex 22 - 28: All other models
DynAdapt	F10				
Highlander	FS3				
Pacifica LP	FS4				
Sierra	FS1				
Pacifica	FS2			24 - 31	
ShockWave	F20	S3	01 - 08	23 - 30	
Highlander MAX	FS3	H5	10 - 12	22 - 31	
Kinterra	ROM	00	01 - 07	22 - 31	
Kinnex	F14	N2	01 - 07	24 - 30	

PART NUMBER EXAMPLE: If ordering an Agilix Foot Module for a 150 lbs. patient with a right-side amputation, moderate activity level, in a size 23 sandal toe, the part number would be: **F15-00-04A23-SR**

STIFFNESS CATEGORY SELECTION

PRODUCT	WEIGHT lbs. / kg.	IMPACT LEVEL		
		LOW	MODERATE	HIGH
Agilix*	110 - 115 / 44 - 52	1	1	2
	116 - 130 / 53 - 59	1	2	3
	131 - 150 / 60 - 68	2	3	4
	151 - 170 / 69 - 77	3	4	5
	171 - 195 / 78 - 88	4	5	6
	196 - 220 / 89 - 100	5	6	7
	221 - 255 / 101 - 116	6	7	8
	256 - 285 / 117 - 130	7	8	9
	286 - 325 / 131 - 147	8	9	
Highlander MAX	326 - 365 / 148 - 166	9		
	286 - 325 / 131 - 147	8	9	10
	326 - 365 / 148 - 166	9	10	11
	366 - 420 / 167 - 190	10		
Kinterra Kinnex	421 - 480 / 191 - 218	11		
	481 - 500 / 219 - 227	12		
	110 - 115 / 44 - 52	1	1	
	116 - 130 / 53 - 59	1	2	
	131 - 150 / 60 - 68	2	3	
	151 - 170 / 69 - 77	3	4	
	171 - 195 / 78 - 88	4	5	
	196 - 220 / 89 - 100	5	6	
221 - 255 / 101 - 116	6	7		
256 - 275 / 117 - 125	7			

PRODUCT	WEIGHT lbs. / kg.	IMPACT LEVEL		
		LOW	MODERATE	HIGH
ShockWave	110 - 115 / 44 - 52	1	1	2
	116 - 130 / 53 - 59	1	2	3
	131 - 150 / 60 - 68	2	3	4
	151 - 170 / 69 - 77	3	4	5
	171 - 195 / 78 - 88	4	5	6
	196 - 220 / 89 - 100	5	6	7
	221 - 255 / 101 - 116	6	7	8
	256 - 285 / 117 - 130	7	8	
	286 - 325 / 131 - 147	8		

IMPACT LEVELS BY ACTIVITY

LOW	MODERATE	HIGH
Walking	Jogging	Running
Biking	Boxing	Competitive Sports
Bowling	Hiking	Football
Camping	Martial Arts	Power Lifting
Dancing	Skating	Snowboarding
Fishing	Skateboarding	Soccer
Gardening	Skiing	Wakeboarding
Golfing	Tennis	
Surfing	Volleyball	

**Agilix:*
 - Not recommended for High Impact.
 - Not recommended for 110 - 115lbs. / 44 - 52kg.

FTC -
 _
_
_
_
_
_
_
_
4 -
 _
_

PRODUCT	STANDARD FOOT SHELLS	GEN_PROFILE - CAP	SIZE	SKIN TONE	TOE	SIDE
Agilix DynAdapt Highlander Highlander MAX Pacifica* Pacifica LP Sierra Kinterra	GENERATION 3: Mid Profile, No Cap	3M - 1	22 - 31*	L: Light M: Medium D: Dark	R: Regular	L: Left R: Right
	GENERATION 2: Sandal Toe, No Cap	2M - 1	22 - 28*		S: Sandal	
ShockWave	GENERATION 3: PRO REGULAR TOE Mid Profile, No Cap	3F - 1	23 - 30		R: Regular	
	GENERATION 2: PRO SANDAL TOE Sandal Toe, No Cap	2F - 1	23 - 28		S: Sandal	
Kinnex	GENERATION 3: Mid Profile, No Cap	3K - 1	24 - 30		R: Regular	
	GENERATION 2: Sandal Toe, No Cap	2K - 1	24 - 28		S: Sandal	

* Pacifica is not available in sizes 22 - 23.

PART NUMBER EXAMPLE: If ordering an Agilix Foot Shell for a 150 lbs. patient with a right-side amputation, moderate activity level, in a size 23 sandal toe, and medium skin tone, the part number would be:
FTC-2M-123M4-SR

AVAILABLE UPON REQUEST

PRODUCT	OTHER FOOT SHELLS	GEN_PROFILE - CAP	SIZE	SKIN TONE	TOE	SIDE
Agilix DynAdapt Highlander Highlander MAX Pacifica* Pacifica LP Sierra Kinterra	GENERATION 3: Mid Profile, With Cap	3M - 0	22 - 31*	L: Light M: Medium D: Dark	R: Regular	L: Left R: Right

* Pacifica is not available in sizes 22 - 23.

REV B 10/2021