





# Heavy Duty Adapter Kit for the Freedom Highlander MAX

	Description	Part Number	Weight Rating	Material
	Heavy Duty Adapter Kit (includes two clamps, one 4-ear connector, one tube)	KIT-00-20200-00	226 kg	Titanium, Graphite
	Tube Clamp 34 mm	ACC-00-20201-00	226 kg	Titanium
	4-Ear Connector for 34 mm System	ACC-00-20202-00	226 kg	Titanium
	Carbon Pylon/Tube	ACC-00-20203-00	226 kg	Graphite