



# MEET THE NEW **KID** ON THE BLOCK

We all know kids can be tough on prosthetic feet; we've got the solution for you! Introducing the ultra low-profile **RUSH Kid**. The all-new **RUSH Kid** features all the same flexibility and durability offered by the RUSH Foot Collection, with a low build height of just 2"! Same smooth roll-through, same dynamic toe-off, most durable pediatric foot available!

#### RUSH Kid PROFILE INFORMATION

<b>Weight Limit</b>	135lbs / 61kg	<b>Heel Rise</b>	3/8" / 10mm
<b>Size</b>	19cm – 21cm	<b>†Suggested Codes</b>	L5981, L5986
<b>Foot Shell</b>	Light or Dark	<b>Warranty</b>	18 months
<b>Product Weight*</b>	505g	<b>Build Height**</b>	2.0" / 51mm



† The listing of this code should not be considered a guarantee for coverage or a payment.  
\* Product weight based on size 20cm, Cat 4 foot module with spectra sock and foot shell.  
\*\*Build height based on size 20cm, Cat 4 foot module with spectra sock, foot shell, and 10mm heel rise.

#### CATEGORY CHART

Weight	lbs	0 - 30	31 - 45	46 - 60	61 - 75	76 - 90	91 - 105	106 - 120	121 - 135
	kg	0 - 14	15 - 20	21 - 27	28 - 34	35 - 41	42 - 48	49 - 54	55 - 61
<b>Low Activity</b>		1P	2P	3P	4P	5P	6P	7P	8P
<b>Moderate Activity</b>		1P	2P	3P	4P	5P	6P	7P	8P
<b>High Activity</b>		1P	2P	3P	4P	5P	6P	7P	8P
<b>High Impact</b>		2P	3P	4P	5P	6P	7P	8P	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

Call us toll-free at 855.450.7300 or visit us at [proteorusa.com](http://proteorusa.com)

f t YouTube in p #HumanFirst



Human First