

MORE PERSONALIZED

PROTEOR
QUATTRO[™]
MICROPROCESSOR KNEE

As the shortest MPK on the market, PROTEOR QUATTRO allows patients to choose from a wider selection of ankles/feet giving them more of what they *need* to live the life they *want*.

	SIZE	22-25cm	26-28cm	29-31cm
Kinterra		13.20"	13.23"	13.38"
ShockWave*		15.44"	15.94"	16.38"
Agilix		14.03"	14.75"	15.30"
Highlander		12.30"	12.84"	13.14"
Pacifica**		11.05"	11.10"	11.25"
Pacifica LP		10.40"	10.64"	10.79"
Sierra		13.81"	14.70"	15.26"
DynAdapt		14.62"	15.05"	15.50"
	SIZE	22-24cm	25-27cm	28-30cm
ROGUE 2/EVAQ8		14.63"	15.13"	15.25"
RAMPAGE/EVAQ8		15.13"	15.63"	15.75"
HiPro/EVAQ8		15.00"	15.38"	15.63"
RAMPAGE LP/EVAQ8		13.00"	13.25"	13.25"
ROVER		10.63"	10.86"	11.00"

**ShockWave 23-30cm / **Pacifica 24-31cm*

Build Heights Based On:

- QUATTRO: 8.5"
- Ankle/Foot: Size 23cm/26cm/29cm, Cat 4 foot module, foot shell, Spectra™ sock, and 10mm heel height

*****Final Build Height Dependent On Socket and Componentry*****



1236 West Southern Avenue #101, Tempe, AZ 85285 855.450.7300

shop.proteorusa.com



#QUATTROKnee
#HumanFirst



Trial QUATTRO and
experience more