

# RUSH™

## RAMPAGE



## Unbridled Capabilities.

The latest addition to the family, the **RUSH Foot RAMPAGE**, has arrived. We've taken the smooth roll-through, increased comfort, and all-terrain capability of its predecessors and created a new, more cosmesis-friendly option, with a lower overall product weight. The **RUSH Foot RAMPAGE** is the perfect solution for individuals seeking a blend of both cosmetic, and high-performance capabilities found only in the RUSH Foot Collection. Are you ready to **RAMPAGE?**

- Cosmesis-ready to meet your aesthetic preferences.
- Lower overall product weight.
- Smooth roll through with unparalleled energy return.
- Biometric ankle motion simulates natural gait.
- Extreme reliability for any use in any environment.
- Manufacturer recommended coding based on AOPA testing guidelines L5987 and L5986.

### RUSH RAMPAGE PROFILE INFORMATION

<b>Weight Limit</b>	365lbs./165kg	<b>Heel Rise</b>	3/8" / 10mm	<b>Size</b>	<b>Build Height**</b>
<b>Size</b>	22cm - 29cm	<b>Warranty</b>	36 Months	<b>22 - 24cm</b>	6.625" / 168mm
<b>Foot Shell</b>	Light or Dark	<b>Suggested Codes</b>	L5987, L5986	<b>25 - 27cm</b>	7.125" / 181mm
<b>Product Weight*</b>	850g			<b>28 - 29cm</b>	7.25" / 184mm

† The listing of this code should not be considered a guarantee for coverage or a payment.  
 \* Product weight based on size 26cm, Cat 4 foot module with spectra sock and foot shell.  
 \*\*Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell, and 10mm heel rise.

### CATEGORY CHART

Weight	lbs	0 - 105	106 - 140	141 - 175	176 - 210	211 - 245	246 - 280	281 - 315	316 - 350	351 - 365
	kg	0 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166
<b>Low Activity</b>		1	2	3	4	5	6	7	8	9
<b>Moderate Activity</b>		1	2	3	4	5	6	7	8	9
<b>High Activity</b>		1	2	3	4	5	6	7	8	9
<b>High Impact</b>		2	3	4	5	6	7	8	9	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

Call us toll-free at 855.450.7300 or visit us at [proteorusa.com](http://proteorusa.com)

f t YouTube in p #HumanFirst



HumanFirst