



Smaller.  
Slimmer.  
Lighter.

# RUSH

## ROGUE<sup>2</sup>

A red circular logo with the text "PDAC" in white at the top and "VERIFIED" in white at the bottom, following the curve of the circle.

The latest innovation in the world's leading fiberglass **RUSH Foot**<sup>®</sup> collection, the **RUSH ROGUE**<sup>®</sup> 2, is here!

With a reduced overall product weight of **10.5%**, **RUSH ROGUE 2** incorporates progressive design features built off the success of the original **RUSH ROGUE** design, delivering **240%** more vertical shock relief at heel strike and **88%** more axial rotation, all in the redesigned vertical shock rotator that's **32%** smaller in volume!

In addition to the traditional foot module, **RUSH ROGUE 2** is also available with the EVAQ8 elevated vacuum system and the H2O water solution, with all three products being PDAC Verified!

### Smaller. Slimmer. Lighter.

The **RUSH ROGUE 2** is unmatched on every level providing an optimal vertical shock experience, axial rotation for improved comfort, predictable performance, ultimate fiberglass strength, and flexibility for daily life and high impact activities.





**Redesigned Vertical Shock Rotator**  
is 32% smaller in volume resulting in a 10.5% reduction in overall product weight & a slimmer profile reducing fatigue during activities.

**Added Axial Rotation**  
reduces socket pressure while improving comfort.

**Increased Vertical Shock Deflection**  
reduces impact forces on the residual limb.

**Solid Keel**  
allows  $\pm 10$  (degrees)  
inversion / eversion

**Heel Shock Bumper**  
provides >90%  
energy return.

**Half Rocker Sole Plate**  
has continuous contact  
with ground for smooth  
rollover, eliminating  
dead spots.

**Long Toe Lever**  
gives better toe response,  
proprioception, &  
energy return.

**Mid-Stance Support**  
provides greater stability &  
a more natural feeling.



**THE ONLY PDAC VERIFIED FIBERGLASS FOOT:**  
L5987 + L5984 + L5781 (RUSH ROGUE® EVAQ8)



## FEATURES

<b>Vertical Shock Rotator</b>		6.9in <sup>3</sup>	10.2in <sup>3</sup>
<b>Vertical Shock Relief</b>		4.25mm	1.25mm
<b>Axial Rotation (each direction)</b>		±15°	±8°
<b>Product Weight*</b>	FOOT	1039g	1155g
	H2O**	883g	1000g
	EVAQ8	1048g	1165g
<b>PDAC Verified</b>	FOOT	L5987 + L5984	L5987 + L5984
	H2O	L5987 + L5984	X
	EVAQ8	L5987 + L5984 + L5781	X
<b>Part Numbers</b>	FOOT	ROG2	ROG***
	H2O	H2R2	H2R***
	EVAQ8	EVQR2	EVQR***

\*Product weight based on a size 26cm, Cat 4 Foot Module with foot shell and Spectra Sock.

\*\*H2O weight and build height does not include Foot Shell or Spectra Sock.

\*\*\*The original RUSH ROGUE®, RUSH ROGUE H2O, and RUSH ROGUE EVAQ8 are no longer available as of July 1, 2021.

## THE ONLY PDAC VERIFIED FIBERGLASS FOOT

