

NON-STERILE FACE MASK

SINGLE USE ONLY



WITH ELASTIC EAR LOOP

PREMIUM QUALITY

PRODUCT DESCRIPTION:

This mask has an integrated nose band with soft rounded earloops and a comfortable inner facial lining.

PRODUCT SPECIFICATIONS:

Protection Level: Level 2 protection

Filtration %: BFE over 98%

Material: 3 layer non-woven polyethelyene fabric

Product SKU: 7000-1689

UOM: 50pcs/bag, 2000pcs/case

Fiberglass Free, Latex Free



SMILESTREAM

SmileStream LLC

www.smilestreamolutions.com

(888) 261-9933

Disclaimer: Pricing and offers subject to change without notice