

UTERINE SOUND & DILATOR



DISPOSABLE UTERINE SOUND

Uterine Sounds (DS) are used for two purposes, one to dilate the cervix and the second to measure the cervix. Sounds are graduated in centimeters and are malleable so they can be shaped to follow the true curvature of the cervix and uterus.

- Latex free, sterile, disposable and ready to use
- Centimeter depth markings
- Sturdy, yet flexible

US-01-025

Box of 25



DISPOSABLE FLEXIBLE DILATOR/OS FINDER

- Length 21.5cm
- OD Size: 2mm - 4mm taper

OS-01-025

Box of 25

Volume discounts and price matching available.
Free samples provided upon request.

1-800-620-2101 www.mpmmedicalsupply.com



265 Willow Brook Road Unit 10, Freehold, NJ 07728 Phone (800) 620-2101 Fax: (800) 620-2127
Email: Sales@mpmmedicalsupply.com


HSG CATHETERS & TRAYS



Your Premium Source
for Quality OBGYN
Products & Instruments



Visit Us Online!
www.mpmmedicalsupply.com




HS CATHETER

The HS Catheter is ideal for Hysterosalpingography and Sonohysterography procedures.

- Latex free, sterile, disposable and ready to use
- Placement sheath decreases the need for a tenaculum
- Patient ambulatory once instrument properly inserted
- Allows for clear view of external cervical Os
- Free dilator with each box of HSG Catheters

TMI1154 5fr Box of 10

TMI1155 7fr Box of 10



SHAPEABLE HS CATHETER

The "Shapeable" HS Catheter is ideal for Hysterosalpingography and Sonohysterography procedures.

- Latex free, sterile, disposable and ready to use
- Placement sheath decreases the need for a tenaculum
- Integrated stylet is malleable and can be re-shaped as necessary
- Angle can be set to align with cervical canal
- Stiffness beneficial for restricted or stenotic cervical canal
- For use with saline or water-based contrast media

TMI1185 5fr Box of 10

TMI1187 7fr Box of 10




SONO-INJECT CATHETER

Ideal for Hysterosalpingography and Sonohysterography procedures.

- Latex free, sterile, disposable and ready to use
- 5fr size for easy injection
- Placement sheath decreases the need for tenaculum
- Allows for clear view of external cervical Os

TMI1170 Box of 10



EZ-INJECT SONOHYSTEROGRAPHY CATHETER

The EZ-Inject SonoHysteroGrapy Catheter can be used for Hysterosalpingography and Sonohysterography procedures.

- Latex free, sterile, disposable and ready to use
- Simple to use – insert and inject
- No distal balloon to block your view of the uterus
- Specially designed top helps prevent leakage
- Cost effective

TMI1162 Box of 10



HSG PROCEDURE TRAY

The Disposable HSG Procedure Tray is ideal for Hysterosalpingography and Sonohysterography procedures.

- Latex free, sterile, disposable and ready to use
- Reduces time spent gathering products
- Ease of ordering, order one product for the procedure
- Easily include HS Catheter with your order

TRAY INCLUDES:

- (1) Wrap, CSR 30"
- (4) Gauze 4" x 4"
- (1) Medium Speculum; Welch Allyn
- (1) Needle, 18G x 1 1/2"
- (2) Swabsticks, OB/GYN 8"
- (1) PVP Ointment, 1 Gram
- (2) PVP Solution, .75 oz
- (1) Disposable Flexible Dilator/OS Finder
- (1) Syringe, 20cc
- (1) Lubricating Jelly
- (1) Drape w/3" Fenestration
- (2) Sponge, Prep Stick 6"
- (1) Panty Liner
- (1) Outer Tray
- (1) Inner Tray
- (1) Extension Tubing Set 36" MF

TMI1157 Tray only Box of 10

TMI1160 with 5 Fr HSG Catheter Box of 10

TMI1161 with 7 Fr HSG Catheter Box of 10