

# SHEERVISION® SV FACE-SHIELD SYSTEM

## SPECIAL INTRODUCTORY OFFER

### ONLY \$49.95 | NOW SHIPPING

In response to the COVID-19 pandemic, SheerVision has designed and launched the SV Face-Shield System. This system is specifically designed to provide face protection for dental / surgical loupe & portable LED headlight wearers.

The innovative arched design of the shield allows the wearing of loupes & portable LED headlights by a surgeon, dentist, dental hygienist, or veterinarian.

The SV Face-Shield System features a transparent, durable, multi-use protective face-shield that includes an adjustable soft headband platform for all-day comfort.

#### **Protective Face-Shield Specifications:**

- Surgical / Dental Loupe & Portable LED Headlight Compatible
- All-Day Comfort through Padded Soft Headband Design
- Protective Premium Full-Face Shield for Medical & Dental Use
- Optically Clear, Transparent, Scratch Resistant, Non-Glare Lens
- Reusable Shield - Thick, Medical Grade PET-G Material
- Recyclable Shield Materials
- Headband Design with Adjustable Knob - One Size Fits All
- Provides Face, Neck, and Ear Coverage
- California Medical Device Manufacturer (CDPH License #95420)
- Reduces risk of contact with dangerous biological materials
- Properly disinfect after each use (70% Alcohol Solution)
- Made in USA
  - Please note loupes & headlights are sold separately.

#### **SV Face-Shield Ordering Information:**

<u>Part #:</u>	<u>Description:</u>	<u>Price:</u>
SV-FS-PREMIUM-SYS	SV-Face-Shield System (Headband w/Shield)	\$ 49.95
SV-FS PREMIUM-4PK	Premium Replacement Shields – 4 Pack	\$ 34.99
SV-FS-PREMIUM-12PK	Premium Replacement Shields – 12 Pack	\$ 79.99
SV-TOP-REPL-4PK	Foam Top Covers – 4 Pack	\$ 5.95
SV-TOP-REPL-12PK	Foam Top Covers – 12 Pack	\$ 16.95
SV-REPL-PAD-2PK	Replacement Headband Pads – 2 Pack	\$ 15.95
MULTI-MINI-CLIP	SheerVision LED Headlight Multi-Clip Mount	\$ 29.95

<https://www.sheervision.com/products/sv-face-shield>



SV Face-Shield Shown w/TTL Loupes & V-RAY LED with Multi-Clip



SV F-S Shown w/Flip-Up Loupes & Ultra-Mini LED on Multi-Clip



SV Face-Shield Shown w/Optional Foam Top Cover, TTL Loupes & Ultra-Mini LED on Multi-Clip

#### SheerVision, Inc.

4030 Palos Verdes Drive N, #104, Rolling Hills Estates, CA 90274 U.S.A.  
 US Toll Free: (877) 678-4274 | International + 1 310 265 8918  
 customerservice@sheervision.com | www.sheervision.com

