

# SHEERVISION® LOUPE FRAME & Rx OPTIONS



**Big Shot**



Black Carbon  
Rx: +2.00 to -2.00  
Single Lens / Bi-Focal Lens (Flip-Up)

**TTL & FLIP-UP**



**Zone XL**



Black Graphite Crimson  
Rx: +2.00 to -2.00  
Single Lens / Bi-Focal Lens

**TTL ONLY**



**Propel**



Black  
Rx: +2.00 to -2.00  
Single Lens / Bi-Focal Lens (Flip-Up)

**TTL & FLIP-UP**



**Clip-On Loupe**

Black  
Rx: Attaches to client's current  
prescription eyeglasses

**CLIP-ON ONLY**



**Performance Titanium - Oval Shape**




Black Platinum Purple (Please contact us for current availability)  
Rx: +12.00 to -12.00  
Single Lens / Bi-Focal Lens

**TTL & FLIP-UP**



**Performance Titanium - Square Shape**



Black Platinum Purple (Please contact us for current availability)  
Rx: +12.00 to -12.00  
Single Lens / Bi-Focal Lens

**TTL & FLIP-UP**

Note: The TTL Rx shown is for the carrier lens. For adding a prescription for near in the back of the barrel, the range is +6.50 to -6.50

**FOR MORE INFO PLEASE ASK YOUR REP OR VISIT SHEERVISION.COM**