

# Home Care Beds

## Accident Prevention

Roscoe Medical's line of Accident Prevention products help protect against injury when used properly. Our floor mats help prevent injury to the patient in the event of accidental roll-off from the bed. Constructed from high-compression foam, the mats provide a firm support surface to be placed beside the bed to create a safer environment for the at-risk patient. Roscoe Bed Rail Pads protect against bumps and bruises that can occur when skin comes in contact with metal bed rails.



30036



70095-BF

### Floor Mats

- Help reduce the possibility of impact-related injuries in the event of an accidental fall from the bed
- Easy-to-clean vinyl cover
- Made of two-inch high density, high compression polyurethane foam
- Encapsulated in a durable Masongard™ cover without zippers (7095 and 7095-BF)
- Flame retardant meets with California Technical Bulletin #117
- Two styles to choose from: One-Piece and Bi-Fold
- Folded sections are held together with hook & look fasteners
- One year warranty



### Full-Length Bed Rail Pads

- Designed for added cushioning and protection
- Deter patient from unassisted bed exit
- Vinyl cover wipes clean for easy disinfecting
- Pads do not interfere with the operation of the bed rails, and fold easily for storage
- Features three hook and loop straps for fast and easy application and removal
- Sold as a pair
- One year warranty

Item #	Description	Color	Size	Qty
30036	Roscoe High Density Foam Floor Mat, Bi-Fold	Blue	24" x 66" x 2"	each
7095	Mason's SafetyCare™ High Density Foam Floor Mat, One Piece	Blue	24" x 66" x 2"	each
7095-BF	Mason's SafetyCare™ High Density Foam Floor Mat, Bi-Fold	Blue	24" x 66" x 2"	each
30038	Roscoe Full-Length Bed Rail Pads	Gray	48" x 15" x 1"	pair



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com