

Upgrade to **Dual Performance**

The new 3M™ Attest™ Dual-reader Biological Indicator System with any-well technology for steam and vaporized hydrogen peroxide

Steam VH2O2

results in 24 minutes





The Power of 2-in-1

Improve your efficiency of processing biological indicators (BIs) with 3M's latest innovation, so you can do more with less. Our 3M™ Attest™ Auto-reader technology incubates and delivers BI results for both steam and vaporized hydrogen peroxide sterilization in one reader.

Any-Well Technology

The new 24-minute Dual System works with Any-Well technology, so both Attest™ BIs 1491 or 1492V (for steam) and 1295 (for vaporized hydrogen peroxide sterilization [VH2O2]) can be incubated in any well of the reader at the same time.

An Upgrade that Gives You More

Get the peace of mind and confidence that comes from knowing every patient receives equipment with a dual Auto-reader system that's

Simple

Current customers can have their existing 3M™ Attest™ Auto-reader model 490 (steam) or 490H (VH2O2) upgraded to start incubating both steam and VH2O2 biological indicators in one unit.



Consistent

Get 24-minute biological indicator results for both steam and vaporized hydrogen peroxide (VH2O2)—providing the ease of consistency.

24 minute BI results*

Familiar & Cost-Sensitive

Use the same trusted 3M[™] Attest[™] Systems with no change to 3M[™] BIs or test packs. The upgrade is no charge, with no new hardware to buy.



Connected

Use the same 3M web application you use now for real-time BI status and results and for digital record keeping.

Your reader will still connect to leading instrument tracking software systems.



To get started with an upgrade of your current 3M™ Attest™ 490 or 490H Auto-reader, contact your 3M Representative, or sign up at go.3M.ca/AttestUpgrade.

Current 3M™ Attest™ Auto-reader 390 users or new 3M customers can contact a 3M sales rep to purchase a new model 490 Auto-reader pre-configured with dual functionality.

Steam Sterilization Monitoring, Now Optimized

3M™ Attest™ Super Rapid Readout Biological Indicators use the same technology you've trusted for years, now optimized to deliver an even faster result.

ISO and FDA Compliant

3M™ Attest™ Super Rapid Readout Bls meet the performance requirements of ISO 11138-1:2017 and ISO 11138-3:2017 and satisfy FDA requirements for biological indicators.

Trusted Steam Sterilization Monitoring Portfolio

- Get the performance of the trusted technology you are familiar with by using 3M™ Attest™ Super Rapid Readout process challenge devices (PCDs).
- The 3M[™] Attest[™] System is the most utilized steam sterilization test pack technology in the U.S., and is now available with the speed of 24-minute results.
- 3M[™] Comply[™] Chemical indicators, Bowie-Dick test packs and indicator tape provide the core four steam sterilization monitoring products.

Steam results in just

minutes*



*When used with 3M™ Attest™ Auto-reader 490 or the 3M™ Attest™ Auto-reader 490H with software version 4.2.7 or greater.

Know. Then Go.

Now every load monitoring is not just possible, it's practical. Take control with 24-minute results and get the peace of mind that comes from knowing that every patient will receive a high standard of care.








Results You Can Count on for Vaporized Hydrogen Peroxide Sterilization

The 3M™ Attest™ System for VH2O2 has the broadest range of FDA clearances in the industry for VH2O2 sterilizers and cycle types.** This provides a truly independent monitor for hydrogen peroxide sterilization assurance needs. The 1295 Biological Indicator system is FDA cleared for the following:



For use with STERRAD® 100S, STERRAD® 100NX (Standard, Flex, Express and Duo cycles), STERRAD® NX (Standard and Advanced cycles), STERRAD® 100NX with ALLClear™ Technology (Standard, Flex, Express and Duo cycles), STERRAD® NX with ALLClear™ Technology (Standard and Advanced cycles).

Amsco® V-PRO® Low Temperature Sterilization System

For use with V-PRO® maX (Lumen, Non Lumen, and Flexible cycles), V-PRO® 1 plus (Lumen and Non Lumen cycles), V-PRO® 1 (Lumen cycles), V-PRO® 60 (Lumen, Non Lumen, and Flexible cycles).

VH2O2 results in just 24 minutes*



3M Provides a Complete VH2O2 System for Monitoring

- Use the key three monitoring products from 3M Bls, Cls, and indicator tape — to leverage the full system of sterilization assurance products.
- A consistent blue-towards-pink colour change for all VH2O2 indicator products, including the chemical indicator (CI) on the biological indicator, internal process CI and external indicator tape, makes it easy to read and interpret results for all of your VH2O2 monitoring needs.

**As of 5/1/2018

The 3M™ Attest™ Brand Advantage



VH2O2



The 3M™ Attest™ Brand advantage gives you a standardized, simplified and supported system for sterilization across steam, VH2O2 and EO modalities.

Fight SSIs from Every Angle

Our 3M portfolio of sterilization assurance products and sterilization monitoring expertise is your first line of defense in your fight to help protect patients from surgical site infections (SSIs).



Load Control Pack Control Exposure Control Equipment Control

Product	Description	Packaging	Shipping
490/490H	3M™ Attest™ Auto-reader 490/490H	1 unit/box	1 box/case
41482V: 41482VF	3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack	24/box	1 box/case
1496V: 1496VF	3M™ Attest™ Super Rapid Readout Steam Challenge Pack	24/box	1 box/case
1491: 1492V	3M™ Attest™ Super Rapid Readout Biological Indicator	50/box	4 boxes/case
1243 (A,B or RE)	3M [™] Comply [™] (SteriGage [™]) Chemical Integrators	500/bag	2 bags/case
		100/bag	10 bags/case
1355	3M [™] Comply [™] Indicator Tapes for disposable wraps	18 mm roll	28 rolls/case
		24 mm roll	20 rolls/case
1322	3M [™] Comply [™] Indicator Tapes	18 mm roll	28 rolls/case
		24 mm roll	20 rolls/case
00135LF	3M [™] Comply [™] Bowie-Dick Plus Test Pack with an early warning sheet	20 packs/box	4 boxes/case
1233	3M™ Comply™ Bowie-Dick Test Pack	6 packs/bag	5 bags/case



Product	Description	Packaging	Shipping
490/490H	3M™ Attest™ Auto-reader 490/490H	1 unit/box	1 box/case
1295	3M™ Attest™ Rapid Readout Biological Indicator	30/bag	4 bags/case
1248	3M [™] Comply [™] Chemical Indicator Strip	250/bag	4 bags/case
1228	3M [™] Comply [™] Indicator Tape	1.9 cm roll	24 rolls/case

To get started with an upgrade of your current 3M™ Attest™ 490 or 490H Auto-reader, contact your 3M Representative, or sign up at go.3M.ca/AttestAdvantage.



Available in Canada from 3M Infection Prevention Division 3M Canada P.O. Box 5757 London, Ontario N6A 4T1 Canada 1-800-563-2921 3M.ca/HealthCare

3M Health Care 2510 Conway Ave St. Paul, MN 55144 USA