Quad-Tray Ultimate

Posterior

The Ultimate in Impression Accuracy

The **Quad-Tray**[®] **Ultimate™ Posterior** is dentistry's most advanced dual-arch metal impression tray. Unlike flexible plastic impression trays, the new Quad-Tray Ultimate Posterior tray will not distort to ensure the ultimate in impression accuracy. Incorporate Quad-Tray Ultimate into your practice and reduce your incidence of impression distortion, costly remakes, and disruption to your office routine.

Ultra-Rigid Aluminum

Fabricated with lightweight, ultra-rigid aluminum, the Quad-Tray Ultimate Posterior is 25% more rigid than its predecessor ensuring distortion-free impressions.

Ultimate Low Profile Distal Bar

At only 2.6mm thick, the low profile distal bar provides the ultimate in prevention of impingement at the retro molar pad. The Quad-Tray Ultimate Posterior tray's distal ends are rounded to provide the ideal fit and comfort.

Ultra-Lock[™] Mesh

At only 0.07mm thick, the Ultra-Lock mesh allows for secure fusion of the preparation and opposing impression, preventing the impression material from separating from the tray during removal.

Recyclable

The Quad-Tray Ultimate Posterior's aluminum construction is 100% recyclable — no more plastic trays ending up in the landfill.



Wide Arch

Accommodating virtually all arch shapes, the wide arch design prevents burn through of the teeth, a common problem with most narrow impression trays.



Low Sidewalls

This design feature prevents palatal impingement and axial roll of the tray wall, a common cause of distortion.





Customizable

Adjust the tray's arch width to suit virtually any patient's arch size.



CRYJT 30 DAY CLARANTE

QUAD-TRAY ULTIMATE POSTERIOR 50-PACK (047501)

QUAD-TRAY ULTIMATE POSTERIOR 250-PACK (047251)



To view the latest news, education and more, visit: www.clinicianschoice.com f in @