


# Exacto Post


FIBRE POSTS

from  angelus<sup>®</sup>



- ✓ Innovative fibres
- ✓ Appropriate polymer matrix
- ✓ Precise adaptation to the canal
- ✓ Highly attractive aesthetic properties

- Root-friendly: elastic modulus similar to dentin;
- High translucency for better aesthetics. Exacto posts have the translucency needed for the favourable aesthetic restorations. The use of innovative fibres in combination with an appropriate polymer matrix confers translucency to the post and highly attractive aesthetic properties to the final restoration;
- Minimally invasive, double taper design;



Highly filled with advanced eGlass fibre for maximum strength

**The shape provides the best of both worlds:**

- A smaller shape inside the canal to maintain as much tooth structure as possible;
- A more robust size in the coronal portion for added strength;
- Extra small size zero (.5 diameter at the tip) for smaller root preparations;
- Excellent radiopacity for better clinical control;
- Pre-silanated and roughened surface for efficient bonding;
- Standardized drill without active tip: Exacto drills exactly reproduce the diameter of the post and do not have an active tip, avoiding the risk of accidental perforations.\*

		Exacto n° 0.5
		Exacto n° 1
		Exacto n° 2
		Exacto n° 3

Drills without active tips have identification colours to match with posts



<b>TRANSLUCENT KIT</b>	856220
#0.5/1/2 - 15 posts	
5x #0.5	
5x #1	
5x #2	
3 drills	
<b>TRANSLUCENT TRIAL KIT</b>	
#0.5 - 5 posts, 1 drill size 0.5	856200
#1 - 5 posts, 1 drill size 1	856201
#2 - 5 posts, 1 drill size 2	856202
#3 - 5 posts, 1 drill size 3	856203
<b>TRANSLUCENT REFILL</b>	
#0.5 - 10 posts	856210
#1 - 10 posts	856211
#2 - 10 posts	856212
#3 - 10 posts	856213

\*Drills come with the kits and the trial kits