

ultrapro™ Tx sweep™



Reduce splatter up to 95%¹

The Ultrapro Tx Sweep disposable prophy angle features a new, innovative brush guard with flexible bristles that prevents saliva and paste from collecting and roping around the cup, which helps prevent splatter. The cup also sweeps the prophy paste back toward the tooth, so the paste is not wasted and lost in the saliva.


Five out of six hygienists experienced a significant reduction in splatter when using the Sweep prophy angle.²

- Brush helps prevent accumulation of saliva and paste on the outside of the cup
- Sweeps paste back toward the tooth, allowing the clinician to continue working and make fewer stops to refill the cup with more prophy paste
- Advanced internal blade design
- Outer ridges for interproximal cleaning
- A 20% shorter head and 25% slimmer neck design gives better access and improved visibility
- Optimal flare
- Comfortable, ergonomic shape
- Not made with rubber latex
- Available in soft or firm cup design


1. Data on file.
2. Data on file.

Splatter Comparison Testing

Sweep prophy angle



Traditional prophy angle



sweep

(Available in soft or firm cup)

GREEN (SOFT) 144-PACK	818357
PURPLE (FIRM) 144-PACK	818358
GREEN (SOFT) 500-PACK	818359
PURPLE (FIRM) 500-PACK	818360