




Our NaviTip™ tips are treated to create flexible tips with a smooth, round end—ensuring they will easily snake through canals to deliver irrigants precisely where they are needed. This treatment can create a golden hue you may see on some of our NaviTip tips.

### NaviTip™ Tips 31 ga with Double Sideport Irrigator

- Double sideports deliver irrigants safely, minimizing the possibility of chemicals being expressed past the apex
- The world's smallest cannula navigates the most intricate canal spaces

Designed for: ChlorCid™, ChlorCid™ V, ChlorCid™ Surf, Ultradent™ EDTA 18% Solution, Consepsis™, and Citric Acid 20%.

755121-21 mm 20pk	755123 - 27 mm 20pk
755122-21 mm 50pk	755124 - 27 mm 50pk



21 mm 27 mm

### NaviTip™ Tips 29 ga and 30 ga

- Provide controlled delivery to the apex
- Flexible, stainless steel cannula easily navigates curved canals
- Available in 17, 21, 25, and 27 mm lengths and 29 or 30 ga cannulae

29 ga delivers paste materials.  
30 ga delivers solutions and gels.



29 ga - 17 mm 29 ga - 21 mm 29 ga - 25 mm 29 ga - 27 mm  
30 ga - 17 mm 30 ga - 21 mm 30 ga - 25 mm 30 ga - 27 mm

	29 ga - 20pk	30 ga - 20pk	29 ga - 50pk	30 ga - 50pk
17 mm - White	755112	751249	751378	751421
21 mm - Yellow	755113	751349	751374	751422
25 mm - Blue	755114	751250	751376	751423
27 mm - Green	755115	751354	751377	751424
Assorted 17-27 mm	755116	751351	751379	-

### Endo-Eze™ Irrigator Tip 27 ga

- Provides ideal reach reducing risk of expressing chemicals past the apex
- Comes with a flexible, blunt cannula with a unique, anti-obturator end
- Non-sterile

Designed for: Ultradent™ 5 ml syringe.

752207 - 25 mm 20pk



### NaviTip™ FX™ Tips 30 ga

- One-of-a-kind brush cleans, scrubs, and irrigates simultaneously
- Rigid 30 ga cannula

Designed for: Ultradent™ Citric Acid 20% Solution.

Listed as an "EXCELLENT" product by a prominent independent research institute.<sup>1</sup>

751452 - 17 mm 20pk  
751454 - 25 mm 20pk

<sup>1</sup>Clinical Research Associates Newsletter, Volume 29, Issue 1, January 2005.



17 mm 25 mm

### Endo-Eze™ Tips

- Work for a variety of endodontic procedures
- Come with flexible, strong cannulae in a variety of sizes
- Bend easily

Designed for: luting materials and air/water delivery.  
Use with: TriAway Adapter, PermaFlo™ DC (20 ga), and other Ultradent syringes.



18 ga 19 ga 20 ga 22 ga

Gauge	Diameter	20pk	100pk
18 ga	0.049"; 1.25 mm	752345	751428
19 ga	0.042"; 1.06 mm	752346	751429
20 ga	0.035"; 0.90 mm	752347	751430
22 ga	0.028"; 0.70 mm	752348	751431

### Capillary Tips

- Evacuates canals and substantially minimize use of paper points
- Narrow, flexible taper accesses curved canals
- Do not use for delivery of irrigating materials or endodontic chemistries

Designed for: Ultradent™ Luer Vacuum Adapter for moisture removal from endodontic canals.

Diameter	20pk	50pk
0.014"; 0.36 mm	752341	753099
0.019"; 0.48 mm	752186	751425

LOK-TITE™



0.014" 0.019"

Ultradent Tips

# DESIGNED TO DELIVER



The chemistries you use are different. Some are chemically activated—needing to be mixed immediately before delivery. Others have varying viscosities. Some work in pits and fissures, some inside canals, and some on enamel surfaces. Each chemistry you use is designed for a specific purpose. Shouldn't the same be true for your tips?

Ultradent makes tips designed to deliver each chemistry we create. Whether you're delivering a flowable composite, a viscous gel, or thick impression material we make the perfect tip for the job. And since our tips are engineered on-site, we test each design to ensure it works perfectly with the chemistry it's intended for.

## Check out our tips with **LOK-TITE™**

Luer Lock tips with Lok-Tite feature double threads that lock the tip into place for increased security—and wings for easy attachment and removal.

### Black Mini™ Tip

- Dispenses large volumes
- Opaque plastic preserves the flow of light-cured materials

Designed for: Ultra-Blend™ plus, Ultradent™ LC Block-Out Resin, PermaFlo™, TriAway Adapter, UltraTemp™, Opalescence™ Boost™, Ultradent™ Diamond Polish Mint, OpalDam™, Opalescence™ Endo, and OraSeal™ Caulking.

751514 - 100pk  
751433 - 500pk

**LOK-TITE™**



### Black Mini™ Brush Tip

- Precise, controlled delivery of aqueous materials
- Tight, adjustable brush fibers minimize bubbles
- Unique to Ultradent

Designed for: Consepsis™, Peak™ SE Primer, Peak™-ZM, Seek™/Sable™ Seek™, Ultradent™ Silane, and Ultradent™ Universal Dentin Sealant.

751169 - 100pk  
751432 - 500pk

**LOK-TITE™**



### Black Micro™ FX™ Tip

- Accommodates various viscosities
- Flocked tip fans out to spread materials in a thin, uniform layer
- 22 ga tip

Designed for: PrimaDry™, and PermaSeal™.

751357 - 100pk  
751434 - 500pk

**LOK-TITE™**



### Black Micro™ Tip · Blue Micro™ Tip

- Pinpoint precision
  - Narrow cannulae accurately deliver materials
  - Black Micro tip: 22 ga, Blue Micro tip: 25 ga
- Black Micro designed for: Ultra-Blend™ plus and Etch Arrest™. Blue Micro designed for: Ultra-Etch™, PrimaDry™, and Opal-Etch™.

751085-Black 100pk    751127 - Blue 100pk  
751435-Black 500pk    751436 - Blue 500pk



### SoftEZ™ Tip

- Tip fibers provide visible, controlled delivery
- Brush fibers facilitate smooth application

Designed for: Enamelast™.

754712 - 50pk  
754711 - 100pk

**LOK-TITE™**



### White Mac™ Tip · White Mini™ Tip

- Dispense large volumes
- All-plastic delivery tips

Designed for: OraSeal™ Caulking, OraSeal™ Putty, Consepsis™ Scrub, and Opalustre™.

751361 - White Mac 100pk  
751480 - White Mac 500pk  
752661 - White Mini 20pk



### Intraoral Tip

- Allows precise placement
  - Attaches to dual-barrel mixing tips
- Designed for: PermaFlo™ DC.

755922 - 20pk



### Ultradent™ Mixing Tip

- Mixes and delivers in one action
- Designed for: UltraTemp™, PermaFlo™ DC, and EndoREZ™.

755920 - 20pk



### Micro 20 ga Tip

- Large-gauge cannula enables consistent flow
- Standard flowable composite delivery tip

Designed for: MTAFlow™, PermaFlo, PermaFlo™ Purple, OpalDam, and OpalDam™ Green.

751252 - 100pk  
751437 - 500pk

**LOK-TITE™**



### SST™ (Surgical Suction Tip)

- Ideal for delicate surgeries
  - Large-diameter tip opening
- Designed for: Ultradent™ Luer Vacuum Adapter for small periodontic or endodontic procedures and controlled suction of Opalescence Boost.

751248 - 20pk

**LOK-TITE™**



### Metal Dento-Infusor™ Tip

- Places hemostatic agents precisely, and effectively removes superficial coagulum
- Blunt, bent 19 ga cannula with padded brush enables gentle pressure on the sulcus
- Ultradent's first tip, the "MDI" remains paramount for successful tissue management

Designed for: Astringedent™, Astringedent™ X, ViscoStat™, ViscoStat™ Clear, and Peak™ Universal Bond.

752559 - 100pk  
752560 - 500pk

**LOK-TITE™**



### Inspiral™ Brush Tip

- Delivers viscous or filled materials smoothly via an internal helical channel and ridge
- Tight, adjustable brush fibers minimize bubbles

Designed for: Composite Wetting Resin, Peak Universal Bond, Ultradent™ Porcelain Etch, Ultra-Etch, UltraSeal™ XT plus, and UltraSeal™ XT hydro.

751123 - 100pk  
751033 - 500pk

**LOK-TITE™**



### Blue Mini™ Dento-Infusor™ Tip

- Offers the same tissue management benefits as the Metal Dento-Infusor tip
- Allows controlled flow of drop-sized quantities
- Outer diameter 1.2 mm

Designed for: Hemostatics and Consepsis.

751086 - 100pk  
751440 - 500pk

**LOK-TITE™**



### Micro Capillary™ Tips

- Bright colors are easily identified against soft tissues
- The world's smallest molded tips (inside diameter: 0.008")

Designed for: Perio materials, endodontics, and Ultradent™ Luer Vacuum Adapter.

751120 - 5 mm 20pk    751121 - 10 mm 20pk  
751802 - 5 mm 100pk    -

**LOK-TITE™**



5 mm    10 mm

