

# When the Jaw Hurts Speed-to-Treat with

**QuickSplint®** 

## Same-day Treatment for Acute Jaw Pain

Jaw and Muscle Sprain/Strain (JAMSS) can lead to chronic orofacial pain, temporomandibular disorders, and headache. It is important for the dentist to diagnose it accurately and treat it promptly – if possible, on the same day. QuickSplint is a same-day immediate splint worn at night to reduce oral habits, reduce forces on the muscles and joints, and encourage healing.

## Diagnose Parafunction & Bruxism

An estimated 85–90% of the general population grind their teeth to a degree at some point during their life. Helping your patient understand and manage parafunction can help prevent tooth wear, tooth fracture, TMJ injury or disorder, associated tension headaches, ear aches, and neck pain. QuickSplint combined with the Parafunction Risk Rating protocol allows you to identify risk and develop a treatment strategy.

## Aid Post-op Healing

Post periodontal surgery, implant placement, third molar extraction, endodontic therapy, or restorative procedure, provide QuickSplint for immediate protection of teeth and restorations and to help the jaw muscles relax and recover.

Highly Rated by Clinicians Report®

An independent, non-profit, dental education and product testing foundation, Clinicians Report®, March 2016



**Make QuickSplint  
Part of Your Practice.  
It's Affordable and Shows  
Your Patients You Care.**

As a Cost-effective Tool, QuickSplint is a smart investment for caring for your patients. QuickSplint uses are reimbursable and conform to dental codes D8210 or D9940 and do not preclude the use of dental code D7880 for a full arch appliance. QuickSplint is FDA 510k cleared for use by a licensed professional.

**Ease of Fabrication  
is a Major Advantage  
of QuickSplint.**

Custom-fit chairside in less than 5 minutes, it is designed to fit a multitude of dental arch configurations and may be fitted on the upper or lower arch as determined by the prescriber.


**Use QuickSplint with any fast-set heavy body VPS.**



**FILL**  
Use fast-set heavy body VPS



**ALIGN**  
on upper or lower arch while patient bites down



**TRIM**  
with scalpel when set. QuickSplint is ready for overnight wear



Same-day treatment shows your patients you care.

Wear marks from clenching and bruxing are recorded on the QuickSplint



**QuickSplint Serves Multiple Purposes**

**General Dentistry**

- Temporary Dental Splint
- Diagnose and Treat Jaw Pain, Tension Headaches, Bruxism
- Diagnose and Treat TMD
- Protect New Dental work
- Minimize Complications
- Implant Placement
- Emergency Splint
- Trial Therapy & Transitional Appliance

**Restorative/Cosmetic Dentistry**

- Case Planning Aid
- Simplify Deprogramming
- Bite Records for Restorations
- Protect Temporaries & New Restorations
- Temporization
- Address Jaw Fatigue
- Transitional Appliance

**Periodontics**

- Easy Periodontal Splint
- Post-Operative Healing Aid
- Implant Placement
- Treat Jaw Pain/Jaw Fatigue

**Endodontics**

- Endodontic Diagnosis
- Preserve Tooth Structure
- Post Treatment Healing Aid
- Address Jaw Fatigue

**Jaw Pain & TMD**

- Treat Acute and Urgent Pain
- Diagnose Parafunction & Bruxism
- Diagnose TMD
- Trial Therapy

**Oral Surgery**

- Treat Jaw Sprain/Strain
- TMD Prevention
- Post-Operative Healing Aid

© 2017 Orofacial Therapeutics, LP.  
QuickSplint® is a registered trademark. All rights reserved.