

# Black Oxide Rubber Dam Clamps

From  


Black Oxide Rubber Dam Clamps are made of special spring tempered stainless steel offering ideal spring action and maximum corrosion resistance.

- These clamps have a lusterless black oxide coating that provides corrosion resistance and reduces eye fatigue by minimizing light reflection
- Produced with rigorous standards and tested for metal fatigue to assure a consistent spring action that will last long-term

## BICUSPID CLAMPS



#0

For small bicuspid and pedo central (#874000)



#2

For lower bicuspid with flat beaks (#874004)



#2A

Similar to #2 but with larger beaks (#874002)



#W2A

Similar to #2 but wingless (#874003)

## MOLAR CLAMPS



#5

For upper molars (#874005)



#56T

For lower molars with festooned beaks (#874017)



#7

For lower molars with flat beaks (#874006)



#8A

For lower molars with strong festooned beaks (#874007)



#W8A

Similar to #8A but wingless (#874008)



#12A

For wisdom teeth with serrated beaks (#874010)



#13A

For wisdom teeth with serrated beaks (#874011)



#14

For half erupted molars (#874012)



#14A

Similar to #14 with larger beaks (#874013)



#26

For upper molars especially for erupted tooth (#874015)

BICUSPID CLAMPS 1 PER PACK	
#0	874000
#00	874001
#2	874004
#2A	874002
#W2A	874003

MOLAR CLAMPS 1 PER PACK	
#5	874005
#56T	874017
#7	874006
#8A	874007
#W8A	874008
#12A	874010
#13A	874011
#14	874012
#14A	874013
#26	874015

LABIAL CLAMPS 1 PER PACK	
#9	874009
#212	874021

## LABIAL CLAMPS



#9

Most popular anterior clamp (#874009)



#212

Cervical clamp for labial cavities (#874021)

# True™ Dental Dam



## Isolate With Confidence.

True™ Dental Dam's ideal surface texture and elasticity allows it to glide smoothly through contacts and tuck into the sulcus for easy placement and optimal marginal seal.

- **Tear Resistant:** Outstanding durability and tear resistance for optimized treatment times.
- **Effortless Placement:** Offered in two pre-cut sizes in medium and heavy gauges to provide ease of placement for any procedure.
- **Isolation and Comfort:** With latex and latex free options, True Dental Dam improves infection control, visibility, and maximizes patient comfort.



Physical Properties	Medium Gauge		Heavy Gauge	
	LATEX	LATEX FREE	LATEX	LATEX FREE
THICKNESS	0.23mm	0.22mm	0.26mm	0.26mm
TENSILE STRENGTH	23.5 MPa	17.5 MPa	23.5 MPa	17.5 MPa
ELONGATION AT BREAK	700%	1100%	700%	1100%
TEAR STRENGTH	85 N/mm	55 N/mm	85 N/mm	55 N/mm



Scan to Read  
Dental Advisor Editor's Choice  
and +++++ Evaluation

<b>LATEX (BLUE) - MEDIUM GAUGE</b>	
5" x 5" 52-pack	838252
6" x 6" 36-pack	838236
<b>LATEX (BLUE) - HEAVY GAUGE</b>	
5" x 5" 52-pack	838253
6" x 6" 36-pack	838237
<b>LATEX FREE (PURPLE) - MEDIUM GAUGE</b>	
5" x 5" 50-pack	838200
6" x 6" 50-pack	838201
<b>LATEX FREE (PURPLE) - HEAVY GAUGE</b>	
5" x 5" 50-pack	838202
6" x 6" 50-pack	838203

“The True Dental dam is easy to place and allows for an excellent seal and isolation.”

- Bob Margeas, DDS

