

True™ Dental Dam



Isolate With Confidence.

True™ Dental Dam's ideal surface texture and elasticity allows it to glide smoothly through contacts and tuck into the sulcus for easy placement and optimal marginal seal.

- **Tear Resistant:** Outstanding durability and tear resistance for optimized treatment times.
- **Effortless Placement:** Offered in two pre-cut sizes in medium and heavy gauges to provide ease of placement for any procedure.
- **Isolation and Comfort:** With latex and latex free options, True Dental Dam improves infection control, visibility, and maximizes patient comfort.



| Physical Properties | Medium Gauge | | Heavy Gauge | |
|---------------------|--------------|------------|-------------|------------|
| | LATEX | LATEX FREE | LATEX | LATEX FREE |
| THICKNESS | 0.23mm | 0.22mm | 0.26mm | 0.26mm |
| TENSILE STRENGTH | 23.5 MPa | 17.5 MPa | 23.5 MPa | 17.5 MPa |
| ELONGATION AT BREAK | 700% | 1100% | 700% | 1100% |
| TEAR STRENGTH | 85 N/mm | 55 N/mm | 85 N/mm | 55 N/mm |



Scan to Read
Dental Advisor Editor's Choice
and +++++ Evaluation

| | |
|---|--------|
| LATEX (BLUE) - MEDIUM GAUGE | |
| 5" x 5" 52-pack | 838252 |
| 6" x 6" 36-pack | 838236 |
| LATEX (BLUE) - HEAVY GAUGE | |
| 5" x 5" 52-pack | 838253 |
| 6" x 6" 36-pack | 838237 |
| LATEX FREE (PURPLE) - MEDIUM GAUGE | |
| 5" x 5" 50-pack | 838200 |
| 6" x 6" 50-pack | 838201 |
| LATEX FREE (PURPLE) - HEAVY GAUGE | |
| 5" x 5" 50-pack | 838202 |
| 6" x 6" 50-pack | 838203 |

“The True Dental dam is easy to place and allows for an excellent seal and isolation.”

- Bob Margeas, DDS

