

Fit Strip™

Flexible
Interproximal finishers with a
Twist



FitStrips are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison Dental Solutions.**

One handy handle! FitStrip's unique attachable handle has been **IMPROVED!** It's even easier to attach (or remove) and makes a world of difference for both clinician ease-of-use and patient comfort. Your hand stays outside the mouth, providing a clear field of view and easier operation.

Do the Twist! Twist the color-coded barrel to adjust FitStrips' curvature and it's automatically locked in place—curved for finishing or straight for IPR. It doesn't get any easier than that!



See our video online!

FitStrip™ makes your procedures easier!

- Orthodontics: The first truly complete IPR/ tooth slenderizing kit
- Composite restorations: Contouring/ finishing all interproximal surfaces
- Crown and bridge: Cement removal and clean-up

Ready for a marathon! FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.

Fit Strip™

Flexible
Interproximal finishers with a
Twist



Your Complete Solution for Interproximal Reduction (IPR)!

NEW! .40mm and .50mm strips easily accommodate maximum reduction cases

NEW! IPR gauges for exact reduction

IMPROVED! Attachable handle for greater ease of use and patient comfort

NEW!

IPR Kit

FPSK05 contains:

- 2 handles
- 1 IPR gauge
- 10 Strips: 2 each .11 mm, 2 each .21 mm, 2 each .30 mm, 2 each .40 mm, 2 each .50 mm

Kits All kits come in a package of 10 FitStrips with 2 handles.

NEW Universal Kit for IPR and finishing/polishing		
1 IPR gauge	FPSK04	\$266
4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)		
6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)		
NEW IPR Kit		
1 IPR gauge	FPSK05	\$266
10 double sided (2 of each .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)		
Starter Kit for finishing/polishing		
2 serrated		
4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)	FPSK01	\$237
4 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)		
Single Sided Kit		
2 serrated		
8 single sided (2 of each .08 mm, .10 mm, .13 mm, .18 mm)	FPSK02	\$237
Double Sided Kit		
2 serrated		
8 double sided (2 of each .11 mm, .15 mm, .21 mm, .30 mm)	FPSK03	\$237

Refills	Single Sided	Double Sided	10 pk w/2 handles	4 ct packs (add - 4 to model #)
Super fine 15 Micron	FPSXFSS 0.08 mm thick	FPSXFDS 0.11 mm thick	Either Refill \$237	\$112
Fine 30 Micron	FPSFSS 0.10 mm thick	FPSFDS 0.15 mm thick	Either Refill \$237	\$112
Medium 46 Micron	FPSMSS 0.13 mm thick	FPSMDS 0.21 mm thick	Either Refill \$237	\$112
Coarse 90 Micron	FPSCSS 0.18 mm thick	FPSCDS 0.30 mm thick	Either Refill \$237	\$112
Medium 46 Micron		FPSTMS 0.40 mm thick	Refill \$237	\$112
Medium 46 Micron		FPSXTMS 0.50 mm thick	Refill \$237	\$112
Serrated	FPSSAW 0.05 mm thick		Refill \$237	\$112
FitStrip replacement handles (2-pack) FPSHANDLE-2			\$41	
IPR Gauge FPSGAUGE			\$55	

*Canadian MSRP

Be sure to check out the website at www.garrisondental.com for info on other innovative Garrison products, demo videos and our 6-month money-back guarantee.

Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456

Offer expires 12/31/2021
To order: Call toll free **888.437.0032**
or online at www.garrisondental.com



Mention Code
INCA121FS

