

**3M™ Imprint™ 4** VPS Impression Material

# Stress less over everyday impressions.

Say "YES" to the VPS impression material with the shortest intra-oral setting time. Move through your day with speed, control and precision.

### Speed.

Reduced chair time and stress for patients—Imprint 4 material has the shortest intra-oral setting time among leading VPS fast-setting impression materials.

#### Control.

Long working time for you. Due to the unique chemistry, the setting time doesn't start until the END of the working time. That means you'll have a sufficient amount of time to load and seat the tray— no stressful race against the clock.

#### Precision.

Hydrophilic. Imprint 4 material is hydrophilic, displacing moisture right from the start in the unset state—which helps prevent detail-disrupting voids.

	Product	Dispensing System	Viscosity  Low High	Setting Version	Maximum Working Time (23°C/73°F)	Maximum Intra-Oral Syringe Time (37°C/98°F)	Intra-Oral Setting Time (37°C/98°F)
Tray Materials	3M <sup>™</sup> Imprint <sup>™</sup> 4 Penta <sup>™</sup> Putty		Putty	Regular Set	1:30	_	2:30
	3M <sup>™</sup> Imprint <sup>™</sup> 4 Penta <sup>™</sup> Heavy 3M <sup>™</sup> Imprint <sup>™</sup> 4 Heavy	<b>► ★</b>	Heavy Body	Regular Set	2:00	_	2:00
	3M <sup>™</sup> Imprint <sup>™</sup> 4 Penta <sup>™</sup> Super Quick Heavy 3M <sup>™</sup> Imprint <sup>™</sup> 4 Super Quick Heavy	<u>►</u> T	Heavy Body	Fast Set	1:15	_	1:15
60	3M <sup>™</sup> Imprint <sup>™</sup> 4 Light	7		Regular Set	_	1:00	2:00
Wash Materials	3M <sup>™</sup> Imprint <sup>™</sup> 4 Super Quick Light	7		Fast Set	_	0:35	1:15
	3M <sup>™</sup> Imprint <sup>™</sup> 4 Regular	7		Regular Set	_	1:00	2:00
~	3M™ Imprint™ 4 Super Quick Regular	7		Fast Set	_	0:35	1:15

## **Ordering Information**

#### Item # Product Information

#### Intro Kits

#### 71494 3M™ Imprint™ 4 Penta™ Heavy Intro Kit

Kit includes: 1 – 300ml base paste; 1 – 60ml catalyst; 1 – 50ml Light cartridge; 1 – 50ml Regular cartridge; 1 – Penta™ cartridge; 10 – Penta™ mixing tips, red; 5 – 3M™ Intra-Oral Syringes, green; 5 – mixing tips, yellow; 5 – intra-oral tips, yellow;

## 1 – VPS tray adhesive 71495 3M™ Imprint™ 4 Penta™ Super Quick Heavy Intro Kit

Kit includes: 1 – 300ml base paste; 1 – 60ml catalyst; 1 – 50ml Super Quick Light cartridge; 1 – 50ml Super Quick Regular cartridge; 1 – Penta™ cartridge; 10 – Penta™ mixing tips, red; 5 – 3M™ Intra-Oral Syringes, green; 5 – mixing tips, yellow; 5 – intra-oral tips, yellow; 1 – VPS tray adhesive

#### 71536 3M™ Imprint™ 4 Penta™ Putty Intro Kit

Kit includes: 1 – 300ml base paste; 1 – 60ml catalyst; 1 – 50ml Super Quick Light cartridge; 1 – 50ml Regular cartridge; 1 – Penta™ cartridge; 10 – Penta™ mixing tips, red; 5 – 3M™ Intra-Oral Syringes, green; 5 – mixing tips, yellow; 5 – intra-oral tips, yellow; 1 – VPS tray adhesive

#### 71534 3M™ Imprint™ 4 Intro Kit

Kit includes: 2 – 50ml Heavy cartridges; 1 – 50ml Light cartridge; 1 – 50ml Regular cartridge; 10 – mixing tips, green; 10 – 3M™ Intra-Oral Syringes, green; 1 – VPS tray adhesive

#### 71535 3M™ Imprint™ 4 Super Quick Intro Kit

Kit includes: 2 – 50ml Super Quick Heavy cartridges; 1 – 50ml Super Quick Light cartridge; 1 – 50ml Super Quick Regular cartridge; 10 – mixing tips, green; 10 – 3M™ Intra-Oral Syringes, green; 1 – VPS tray adhesive

#### 3M™ Penta™ Refills—Tray Materials

Each refill includes: 2 – 300ml base pastes; 2 – 60ml catalysts

71484 3M™ Imprint™ 4 Penta™ Heavy

71485 3M™ Imprint™ 4 Penta™ Super Quick Heavy

71486 3M™ Imprint™ 4 Penta™ Putty

Refills—Tray Materials Each refill includes: 4 – 50ml cartridges; 10 – mixing tips, green

71492 3M™Imprint™ 4 Heavy

71493 3M™ Imprint™ 4 Super Quick Heavy

#### Item # Product Information

Economy Packs—Tray Materials
Each Economy Pack includes: 10 – 50ml
cartridges

71531 3M™ Imprint™ 4 Heavy Economy Pack

71532 3M<sup>™</sup> Imprint<sup>™</sup> 4 Super Quick Heavy Economy Pack

#### Refills—Wash Materials

Each refill includes: 4 – 50ml cartridges; 5 – mixing tips, yellow; 5 – 3M<sup>™</sup> Intra-Oral Syringes, green

71488 3M™ Imprint™ 4 Light

71487 3M™ Imprint™ 4 Regular

71490 3M<sup>™</sup> Imprint<sup>™</sup> 4 Super Quick Light

71489 3M™ Imprint™ 4 Super Quick Regular

#### for 2M™ Por

77784 3M™ Imprint™ 4 Heavy

TTTOT OUT IMPINIT THEAVY

77785 3M™ Imprint™ 4 Super Quick Heavy

77786 3M™ Imprint™ 4 Penta™ Putty

#### Cartridge

or 3M™ Pentamix™ 2 Automatic Mixing Uni

76994 Universal 3M™ Penta™ Cartridge

Accessories (includes: 50 - tips or syringes)

77949 3M™ Penta™ Mixing Tips - Red

71452 3M™ Garant™ VPS Mixing Tips - Yellow

71462 3M™ Garant™ VPS Intra-Oral Mixing Tips – Yellow

71450 3M™ Garant™ VPS Mixing Tips - Green

71506 3M™ Intra-Oral Syringe – Green

Warranty: 3M warrants this product will be free from defects in material and manufacture at the time of purchase. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M's sole obligation shall be repair or replacement of the 3M product. Limitation of Liability: Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.



3M Oral Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249 Web 3M.com/dental



# Customer Care Center: 1-800-634-2249 3M.com/greatimpressions

3M, Garant, Imprint, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH. Please recycle. Printed in U.S.A. © 3M 2021. All rights reserved. 70-2013-0580-5 (Rev. 8/2021)