

PEAK™

UNIVERSAL ADHESIVE SYSTEM



CHOOSE YOUR PATH TO A STRONGER BOND

ETCH AND RINSE



NO RINSE



ULTRA-ETCH™ ETCHANT



PEAK™ SE PRIMER

PEAK™ UNIVERSAL BOND ADHESIVE

AVAILABLE
IN BOTTLES &
SYRINGES!



HIGHEST BOND STRENGTH TO DENTIN AND ENAMEL!¹⁻²

Peak Universal Bond adhesive ensures a strong, lasting bond whether you choose a self-etch or a total-etch technique.

- Bonds to dentin, enamel, porcelain, metal, composite, and zirconia
- Ideal for indirect and direct bonding, as well as post and core procedures
- Works with self-etch and total-etch techniques

All paths lead to
Peak Universal Bond
adhesive.

Peak™ Universal Bond

Light-Cured Resin Adhesive

REFRIGERATE



4554 - Peak Universal Bond Self-Etch Syringe Intro Kit

- 1 x 1.2 ml Peak Universal Bond syringe
- 1 x 1.0 ml Peak SE Primer syringe
- 20 x Black Mini Brush tips
- 20 x Inspiral Brush tips

REFRIGERATE



4551 - Peak Universal Bond Total-Etch Syringe Intro Kit

- 1 x 1.2 ml Peak Universal Bond syringe
- 1 x 1.2 ml Ultra-Etch syringe
- 20 x Blue Micro tips
- 20 x Inspiral Brush tips

REFRIGERATE



4541 - Peak Universal Bond Self-Etch Bottle Kit

- 1 x 4 ml Peak Universal Bond bottle
- 4 x 1.0 ml Peak SE Primer syringes
- 40 x Black Mini Brush tips
- 50 x Mixing Wells
- 50 x Micro Applicator brushes

REFRIGERATE



4542 - Peak Universal Bond Total-Etch Bottle Kit

- 1 x 4 ml Peak Universal Bond bottle
- 4 x 1.2 ml Ultra-Etch syringes
- 40 x Blue Micro tips
- 50 x Mixing wells
- 50 x Micro Applicator brushes

REFRIGERATE



5135 - Peak SE Primer Syringe 4pk

1.0 ml syringes

REFRIGERATE




4553 - Peak Universal Bond Syringe 4pk

4552 - Peak Universal Bond Syringe 20pk

1.2 ml syringes

SUPERIOR BOND STRENGTHS!^{1,2}

■ Dentin Bond Strength (MPa) ■ Enamel Bond Strength (MPa)



1. Peak Universal Bond was reviewed favorably by an independent non-profit dental education and product testing institute, *Clinicians Report*, August 2012. Full data is available upon request. 2. Internal testing data.

*Trademark of a company other than Ultradent Products, Inc.

[†]Total-etch technique [‡]Self-etch technique

REFRIGERATE

4543 - Peak Universal Bond Bottle 1pk

4 ml bottle



4545 - Mixing Wells 100pk



4548 - Micro Applicator Brush 400pk

200 x Each color

