

TECHNICAL DATA SHEET
Dental Gauze
Intended Use

Used in Healthcare Facilities for management of topical wounds and temporary control of external surface bleeding such as minor cuts and minor laceration.

Characteristic	Specification	
Material Composition	Carboxymethyl cellulose	
Description	Soft, glossy, slightly stretchable woven material with cut and/or weaved edges.	
Color	White or Light yellow	
Dimension	Width (mm)	Length (mm)
	19.0	19.0

Additional Information

Storage Condition	Country of Origin
Store in a cool, dry place	China

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment	Marketed in USA as	Marketed in Canada as
4900-STAT	Medicom SafeGauze Topical Hemostatic Dressing	1 Dressing/Blister, 10 Blisters/Blister Pack, 2 Blister Packs/Box, 200 Boxes/Case, 4000 Dressings/Case	1,9 x 1,9 cm (0.75" x 0.75")	White or Light Yellow	4 Ply	Sterile	Unclassified	Not Marketed in Canada
4900-STATC	Medicom SafeGauze Topical Hemostatic Dressing	1 Dressing/Blister, 10 Blisters/Blister Pack, 2 Blister Packs/Box, 200 Boxes/Case, 4000 Dressings/Case	1,9 x 1,9 cm (0.75" x 0.75")	White or Light Yellow	4 Ply	Sterile	Not Marketed In USA	Medical Device Class 2