

**TECHNICAL DATA SHEET**
**NITRILE, NON-STERILE MEDICAL EXAMINATION GLOVE: ASTM D 6319-19**
**Intended Use**

Non-sterile, single use, medical examination gloves are intended to prevent transmission of a wide variety of diseases to both patients and health care personnel. Specifically, gloves are intended to act as a barrier between the healthcare provider and patients to prevent contamination between patient and caregiver.

Characteristic		Acceptable Limit	
Dimension & Tolerance	Width (mm)	Extra-Small (XS)	60 - 85
		Small (S)	70 - 90
		Unisize (S-M)	75 - 95
		Medium (M)	85 - 105
		Large (L)	100 - 120
		X-Large (XL)	110 - 130
		XX-Large (XXL)	120 - 140
Length (mm)		XS & S: $\geq 220$ Others: $\geq 230$	
Thickness (mm)	Palm	$\geq 0,05$	
	Finger		
Physical Requirements	Before ageing	Tensile Strength	$\geq 14$ MPa
		Ultimate Elongation	$\geq 500\%$
	After ageing	Tensile Strength	$\geq 14$ MPa
		Ultimate Elongation	$\geq 400\%$
Powder Residue	Powder-Free Gloves	$\leq 2,0$ mg/glove	
	Powdered Gloves	$\leq 10,0$ mg/dm <sup>2</sup>	
Description	Ambidextrous.		
Color	Refer to Appendix 1.		

**Additional Information**

Storage Condition	Country of Origin
Store under cool, dry place. Avoid excessive heat (Maximum 40°C, 104°F). Recommended storage temperature: 10-30°C. Shield open box from direct sun and fluorescent light.	Malaysia or China or Thailand or Indonesia

**Appendix 1 List of Item Codes and Attributes**

## Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment	Marketed in USA as	Marketed in Canada as
1185	Medicom SafeBasics "True Fit Thin" Powder-Free Nitrile Medical Examination Gloves	300 Gloves/Box, 8 Boxes/Case, 2400 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Country of origin: Malaysia	Non-Sterile	Medical Device Class 1	Medical Device Class 2
1185EX	Medicom SafeBasics "True Fit Thin" Powder-Free Nitrile Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Country of origin: Malaysia	Non-Sterile	Not Marketed in USA	Not Marketed in Canada
1186	Medicom SafeBasics "True Fit Thin" Powder-Free Nitrile Medical Examination Gloves	300 Gloves/Box, 8 Boxes/Case, 2400 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large	Pink	Country of origin: Malaysia	Non-Sterile	Not Marketed in USA	Medical Device Class 2