



**Ready to
Grind!**

A Proven Universal
Composite that is
Perfect for Posterior Teeth.

KURARAY
DENTAL



Features & Benefits

Accurate color matching

High strength and durability

High polishability

Easy handling

Does not stick to instruments

Sufficient radiopacity

Superior physical properties

Available in syringes and PLTs

Very economical per use application

CLEARFIL AP-X RESTORATIVE SYSTEM **CLEARFIL AP-X PLT**

Ten years of clinical experience and scientific research

guarantee **Clearfil AP-X** offers superior aesthetic appearance and exceptional physical properties. Compared to all other available composites, **Clearfil AP-X** series offers 5 good reasons to change.

CLEARFIL AP-X RESTORATIVE SYSTEM CLEARFIL AP-X PLT

1 Low wear and high fracture resistance

Clearfil AP-X series is a micro-hybrid composite consisting mainly of high density fine barium glass particles. The fine barium glass particles make up 85% of the composite's weight. This is one of the reasons for the unusually low wear and high fracture resistance of the **Clearfil AP-X** series; a truly convincing value for durable posterior restorations.

2 Excellent Polishability

The high density fine filler particles (smallest 0.04 um) provides excellent polishability. The aesthetically pleasing result makes **Clearfil AP-X** series the perfect choice for challenging anterior restorations.

3 Low Polymerization Shrinkage

The high quality of our inorganic filler is one of the reasons for a very low polymerization shrinkage of 1.9%; reducing stress forces against cavity walls dramatically. As a result, the risks of post-operative sensitivities and secondary cavities due to leakage are significantly limited.

4 Visible Radiopacity

The composite's high percentage of barium glass particles offers excellent radiopacity for x-ray evaluations.

5 Easy Handling

Clearfil AP-X series is made with a special non-sticky formula, preventing the composite from sticking to any instruments or contact matrix. The material offers convenient adaptation for use with additional composite layers making it an excellent choice for anterior restorations. In addition, the consistency of **Clearfil AP-X** series allows for precise shaping of posterior restorations. Convenient PLT compules improve the application process and limits the possibility of cross-infection.

COMPOSITE RESIN FILLING MATERIAL COMPARISON

Product	Manufacturer	Filler Content (wt%)	Flexural Strength (MPa)	Vickers Hardness	Tooth Brush Wear (mm ³)
AP-X	Kuraray Medical	85.5	208	134	0.20
TPH Spectrum	Dentsply	77.0	131	61	0.92
Tetric Ceram	Ivoclar Vivadent	76.0	164	N/A	N/A
Charisma	Kulzer	69.5	111	61	0.92
Solitaire	Kulzer	65.0	71	56	N/A
Point 4	Kerr	76.0	142	61	0.36
Z100	3M-Espe	84.5	174	119	0.21
Filtek Z250	3M-Espe	77.7	165	84	N/A
Dytract AP	Dentsply	73.0	104	51	1.34
Esthet-X	Dentsply	81.4	136	70	0.53
F2000	3M-Espe	84.0	67	130	0.98
Filtek Supreme	3M-Espe	73.6	126	92	N/A

Tests performed by Kuraray Medical Inc.

CLEARFIL™ AP-X series components

SYRINGES

- #1721-KA CLEARFIL™ AP-X (A2)
- #1722-KA CLEARFIL™ AP-X (A3)
- #1723-KA CLEARFIL™ AP-X (A3.5)
- #1724-KA CLEARFIL™ AP-X (A4)
- #1725-KA CLEARFIL™ AP-X (B2)
- #1726-KA CLEARFIL™ AP-X (B3)
- #1727-KA CLEARFIL™ AP-X (B4)
- #1728-KA CLEARFIL™ AP-X (C2)
- #1729-KA CLEARFIL™ AP-X (C3)
- #1730-KA CLEARFIL™ AP-X (C4)
- #1731-KA CLEARFIL™ AP-X (XL)
- #1732-KA CLEARFIL™ AP-X (C1)
- #1733-KA CLEARFIL™ AP-X (CL)
- #1734-KA CLEARFIL™ AP-X (HO)

PLTs

- #1715-KA CLEARFIL AP-X PLT Intro and Gun kit (A2, A3, B2, XL: 10 tips ea./Clearfil SE BOND*/Clearfil AP-X PLT Gun)
 - #1740-KA CLEARFIL™ AP-X PLT Intro Kit (A2, A3, B2, XL: 10 tips ea./Clearfil SE BOND*)
 - #1741-KA CLEARFIL™ AP-X PLT 20 tips (A2)
 - #1742-KA CLEARFIL™ AP-X PLT 20 tips (A3)
 - #1743-KA CLEARFIL™ AP-X PLT 20 tips (A3.5)
 - #1744-KA CLEARFIL™ AP-X PLT 20 tips (A4)
 - #1745-KA CLEARFIL™ AP-X PLT 20 tips (B2)
 - #1746-KA CLEARFIL™ AP-X PLT 20 tips (B3)
 - #1751-KA CLEARFIL™ AP-X PLT 20 tips (XL)
 - #1754-KA CLEARFIL™ AP-X PLT 20 tips (HO)
- *CLEARFIL™ SE BOND Primer 1.5ml / CLEARFIL™ SE BOND Bond 1.25ml / Outer Case, Accessories

Clearfil AP-X and SE Bond



Clearfil AP-X, is the ideal direct material to use with **SE Bond**.

Clearfil SE Bond

Clearfil SE Bond is a self-etching primer and bonding system that is highly recommended and scientifically proven.

Clearfil SE Bond's water-based primer significantly reduces procedural time while providing excellent bond strength and sealing performance. The unique delivery system inhibits direct finger contact preventing any risk of contact allergies.

Clinical Indications

Class I, II, V restoration of posterior teeth
Class III, IV, V restoration of anterior teeth
Cervical cavities or defects involving root surface

Distributed by:

Kuraray America, Inc.

101 East 52nd Street, 26th Floor
New York, NY 10022

(800) 879-1676

www.kuraraydental.com

Manufactured by:

Kuraray Medical, Inc.

1621 Sakazu, Kurashiki,
Okayama 710-8622 Japan

Kuraray Dental...over 25 years of committed, focused, scientific work and leadership with relentless search for better solutions—products that enhance your skill and care.

KURARAY
DENTAL

World's Leading Manufacturer of Dental Products