

Rally™

Mini Polishers for Composite

Rally gets you to the finish line faster!

The composite finish line that is. With zero reprocessing and sterilization hassles and a lustrous diamond final gloss, Rally mini-polishers get you across the finish line and on to your next patient.

Get to the point – Every time!

Single patient Rally means that you get a clean, fresh, properly shaped polisher every time.

Mini, but mighty!

Rally's small size and tapered shaft create a truly agile polisher capable of reaching the smallest anatomical features without blocking your visibility.

A clean get-a-way!

You don't want to leave a trail of crumbs or spattered polishing paste on the track. Rally's diamond impregnated high shine final polishers won't crumble during use and don't require paste to create a beautiful, durable shine.



Give your mechanic the day off! (Maybe your assistant too!)

Rally is fully assembled – no little mandrels to lose track of, no discs to snap on, there's just nothing you need to mess with – including cleaning and sterilization.




Rally™


Mini Polishers for Composite



**Remove -
Choose the Blue**
Polishers for fast removal
and contouring



**Smooth -
Choose the Red**
Polishers for smoothing
and silky shine



**Shine -
Choose the Gray**
Polishers for ultra
high shine

Item	Description	Qty.	Price
FPCK01	Composite Polisher Complete Kit	60 (10 each)	\$200

Item	Refills	Qty.	Price
FPCC010	Blue coarse grit point composite polisher	30	\$76
FPCC040	Blue coarse grit cup composite polisher	30	\$76
FPCM010	Red medium grit point composite polisher	30	\$76
FPCM040	Red medium grit cup composite polisher	30	\$76
FPCF010	Gray fine grit point composite polisher	15	\$76
FPCF040	Gray fine grit cup composite polisher	15	\$76



Slow speed right angle latch shank. Use with water.



Do not exceed
maximum speed



Wear breathing
protection



Wear eye
protection



Use waterspray
minimum 50ml/min

*Canadian MSRP



Be sure to check out the website at www.garrisondental.com for info on other innovative Garrison products, demo videos and our 6-month money-back guarantee.

Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456

Offer expires 12/31/2021

To order: Call toll free **888.437.0032**
or online at www.garrisondental.com



Mention Code
INCA121RP

