

● ABOUT US ●

Here at Dr Physio USA, we make the most imaginative, unwinding, and restorative, wellbeing and individual care Products.

In the hubbub of today's lifestyle, it becomes impossible to nurture yourself at times. In the chaos of maintaining every aspect of life efficiently, we forget that our bodies are not machines. They need absolute care and tenderness. Pain in neck, shoulder, back or foot should never be taken for granted. They can lead to serious ailments in the long run. To help you ease the situations of painfulness, Dr Physio USA brings to you excellent ways of relaxation and pain relief through our premium gadgets. We offer a non-intrusive answer to help renew and reestablish your muscles and relieve throbbing pain with our various products

The variety of items and massagers we offer are intended to target particular areas of the body viz neck, shoulder, back, foot and so on. the adaptable capabilities of our gadgets are meant for utmost body relaxation after a long day at work. Since we are continually in pursuit of excellence, our research and development team at NY, USA works diligently to ensure that each of our innovations are available at a value you would appreciate.

Sit Back, Relax and enjoy the amazing Dr Physio experience by Dr Trust!!

Dr Physio from Dr Trust
www.drphys.io | www.drtrust.in

**Dr Physio
Head & Scalp Massager-1028**

● INTRODUCTION ●

Dr Physio Head & Scalp Massager-1028, the original 3D meridian massage cluster comes with 6 kinds of massage effect including opening, closing, thrusting, drawing, rubbing, and scratching. It can promote blood circulation, improve scalp troubles pass through the meridians, eliminate fatigue, eliminate dandruff, stop itching, improve sleep. It is designed to provide the most comfortable meridian massage experience.

● SCOPE OF APPLICATION ●

For hair loss in the head, less hair, poor hair quality, the massager has a significant effect and can significantly improve fatigue.

● INTENDED USER ●


1. People who often feel tired.
2. College students, middle school students and primary school students.
5. People who often have headaches and dizziness.
6. People who suffer from insomnia and dreaminess.

USE PROHIBITED

People who have the following conditions should not use:

1. Patients with hemorrhagic head.
2. Patient undergoing head surgery within three months.

● PRODUCT DESCRIPTION ●



● METHOD OF OPERATION ●

1. Press the button to turn it on first, the default is slow, the second is fast, the third is reversed quickly.
2. Each press is repeated for the first time.
3. Press and hold for 3 seconds to shut down.

● CHARGING ●

If red colour's light start flashing, it means the machine need to be charged immediately.

1. During the charging process, the battery indicator light is always red. When the battery is fully charged, the indicator light turns to yellow.
2. Do not use this product during charging (the charging is not supported in the power-on state. The product automatically stops working during charging for self-protection).
3. The charging time is 2 hours.

● POWER ON, MODE SELECTION, AND SHUTDOWN ●

Power On

1. After charging, the massager will start functioning with a short press of power key.
2. The default slow speed starts when powering it ON.
3. The default time is 10 minutes.

Mode Selection

After turning the machine ON, press the switch key briefly to change and choose the massage mode and knead slowly.

- Fast kneading in the forward direction;
- Slow kneading in forward direction;
- Slow kneading in reverse direction;
- Fast kneading in the reverse direction.

● SHUTDOWN ●

1. The massager turned off with a long press of the power switch.

● PRECAUTIONS ●

To ensure that you use the scalp massager correctly and safely please read the instructions carefully.

- Do not use the massager at high temperatures.
- If the head has redness and inflammation; bruises and other discomfort symptoms, or if you have had a head surgery, do not use it.
- Do not use if you have bleeding disorders such as hemorrhagic disease or hemophilia.
- For skin problems and obstacles caused by silicone products, please consult a physician before purchasing.
- Patients with severe stroke hypertension and hypotension should use the massager with caution. Please consult a professional physician. If there is any discomfort during use, stop using it immediately.
- Please avoid using the massager after drinking or when you are not feeling well.
- Please use the power adapter configured with this product and pay attention to electricity safety.
- The massager is designed as a home electronic health care product and is not suitable for diagnostic and therapeutic purposes.
- Do not place the charging base in a humid environment.
- Do not charge in excessively humid conditions.

● COMMON FAULTS AND TROUBLESHOOTING ●

Fault conditions	Cause Analysis	Processing method
Power dead	Is the battery exhausted	Please charge the battery
Massage strength	Is the battery exhausted	Please charge the battery

● CUSTOMER SUPPORT ●

CONTACT ADDRESS

USA
Nureca INC.USA
276 5th Avenue, Suite 704-397,
New York (NY) - 10001, USA
INDIA
Corporate Office (Mumbai)
Nureca Limited
128 Gala Number Udyog Bhavan,
1st Floor Sonawala Lane, Goregaon East
Mumbai City Maharashtra 400063
Contact us
India: +91-7527013265 /+91-9356658436
Website: www.drtrust.in
Corp Website: www.nureca.com
Email: customer@nureca.com
Connect with us on social networks
Facebook: @drtrust
Instagram: @drtrustin
Youtube: NurecaUsa



Scan to View
Product Demo Video
www.drtrustusa.com/1028