

LIPOSUCTION

PLASTIC SURGERY INSTRUMENTS





FLUID INFUSION

Infusion needles

High grade stainless steel. 6 spaced holes in a spiral around the tip. For fluid infusion using leur lock syringes, handles or giving set tubing.



Product Code	Diameter (mm)	Length (cm)
GI-INF/N 215	2	15
GI-INF/N 223	2	23
GI-INF/N 230	2	30
GI-INF/N 315	3	15
GI-INF/N 323	3	23
GI-INF/N 330	3	30

Curved

Product CodeDiameter (mm)Length (cm)GI-INF/N1.623CRV1.623

Infusion handle non valve

Comfortable handle to link tubing to needle, lack of valve reduces leakage. Single shaft allows infiltration and suction.







Infusion handle with Valve

This infusion handle, with on/off button, offers maximum control for infusion during the tumescent technique.









MICRO INFILTRATION CANNULAE

High grade stainless steel for re-injection using Luer lock syringes.

Product Code	Diameter (G)	Length (cm)
GI-MIC16/A7	16G	7
GI-MIC16/A9	16G	9
GI-MIC16/B7	16G	7
GI-MIC16/B9	16G	9
GI-MIC16/B12	16G	12
GI-MIC16/C7	16G	7
GI-MIC16/C9	16G	9
GI-MIC18/A9	18G	9
GI-MIC18/B9	18G	9





Curved Cannula

Product Code	Diameter (G)	Length (cm)
GI-CURV/18B7	18G	7
GI-CURV/16B12	18G	12
GI-CURV/16B15	18G	15

Ultra fine

Ultra fine facial needle with unblocking Stylet

Product Code	Diameter (G)	Length (cm)
GI-235FN/LL	23G (0.7mm)	5
GI-217FN/LL	21G (0.8mm)	7
GI-2185FN/LL	21G (0.8mm)	8.5
GI-15FN/LL	20G (1.0mm)	5
GI-207FN/LL	20G (1.0mm)	7
GI-189FN/LL	18G (1.2mm)	9



Micro Drop

Curved Obtuse tip for accurate placement of micro droplets.

Product Code	Diameter (mm)	Length (cm)
GI-1285MD	1.2	8.5





we can produce excellent products quickly so... Why not design your own unique cannula



DISSECTORS

V-Dissector

Undermining fibrous tissue, widens at the tip and ideal for breast surgery.

Product Code Diameter (G) Length (cm) GI-VDISS167 16G 7 GI-VDISS169 16G 9

W-Dissector

For undermining small areas of fibrous tissue, less aggressive and finer than V Dissector.

Product Code	Diameter (G)	Length (cm)

16G

HAY Tunnel Dissector

For blunt dissection tunneling when a V-dissector may be considered aggressive i.e. around Implants. Slight upward curve to tip in order to inhibit direct forceful contact, comfortable handle with orientation marking.

Product Code	Diameter (mm)	Length (cm)	
GI-HTD	2	12	

GI-WDISS

ADAPTORS

Toomey to luer lock adaptor

Transfer aspirate from 60cc Toomey to any size

luer lock syringe..

Product Code
GI-MIC/LL/S



Luer lock to luer lock adaptor

Transfer between any size luer lock syringes. Product Code

GI-MIC/LL/LL



8.5

Toomey syringe to luer lock needle adaptor

Product Code GI-MIC/LL/N





Open harvest

GI-322MLL

Guarded fully open tip, designed specifically for tissue harvest.

MICRO

HARVEST Tissue harvest using luer lock syringes.

All Cloverleaf luer lock bases have a thumb platform for comfort and orientation.

Product Code	Diameter (mm)	Length (cm)
GI-MIC315/HARV	3	15
GI-MIC322/HARV	3	22



22


30 Hole harvest cannula

Harvest fat in micro droplets for fine graft re-injection

Product Code	Diameter (mm)	Length (cm)
GI-2.515/30HOLE	2.5	15
GI-315/30HOLE	3	15

Submental

One hole, flat backed for areas closer to the skin



Mercedes tip

3 holes spaced evenly around the tip. The most popular tip style.



CELL SELECT CANNULAE RANGE

Level 1 - Gentle Rasp

Product Code	Diameter (mm)	Length (cm)
GI-2515G/21 HOLES	2.5	15
GI-320G/30HOLE	3	20

Level 2 - Medium Rasp

Product Code	Diameter (mm)	Length (cm)
GI-2515G/21 HOLES	2.5	15
GI-320G/30HOLE	3	20

Level 3 - Aggressive Rasp

Autoclavable steel, holds 60 10cc syringes or 12 60cc syringes

Product Code	Diameter (mm)	Length (cm)
GI-2515G/21 HOLES	2.5	15
GI-320G/30HOLE	3	20

The ingenious Cell Select micro harvest cannula range is specifically designed to capture smaller aliquots of healthy fat cells via numerous 1mm micro-ports. The harvest of smaller cell packages results in easier, smoother re-injection through finer injection cannula. The spacing and increased surface area of the apertures per cannula tip results in an effective and constant collection speed, and greatly reduces blocking

3









POWERED SUCTION

One piece cannulae

High grade Stainless Steel shaft, for aspiration with powered suction.





Diameter (mm)	Length (cm)
2	10
3	15
3	25
4	30
5	30
	2 3 3 4

Product Code	Diameter (mm)	Length (cm)
GI-325/30 HolesM	3	25

Product Code	Diameter (mm)	Length (cm)
GI-325LM	3	25
GI-430LM	4	30
GI-530LM	5	30
	•	00

Product Code	Diameter (mm)	Length (cm)
GI-210TSM	2	10
GI-320TSM	3	20
GI-430TSM	4	30

Product Code	Diameter (mm)	Length (cm)
GI-210SM	2	10
GI-315SM	3	15

Product Code	Diameter (mm)	Length (cm)
GI-215OHM	2	15
GI-320OHM	3	20
GI-430OHM	4	30

Mercedes

3 holes spaced evenly around the tip. The most popular tip style.



30 Hole harvest cannula

Harvest fat in smaller parcels, reduces blocking.



Lantern

Chinese lantern shaped tip, useful for fibrous areas and popular in Gynecomastia procedures.



The Scream

Two holes above a single hole, a cannula lookalike of the famous Munch painting.



Submental

One hole, flat backed for areas closer to the skin.



One hole

A single downward facing aperture.





TOOMEY SYRINGE SUCTION



Toomey Cannulae High grade Stainless steel shaft, for aspiration with 60cc Toomey syringes.



Mercedes tip

3 holes spaced evenly around the tip. The most popular tip style.



30 Hole harvest cannula

Harvest fat in smaller parcels, reduces blocking.



Product Code	Diameter (mm)	<u>Length (cm)</u>
GI-210MTS	2	10
GI-315MTS	3	15
GI-318MTS	3	18
GI-325MTS	3	25
GI-430MTS	4	30

Product Code	Diameter (mm)	Length (cm)
GI-2.515/30 HolesTS	2.5	15
GI-320/30 HolesTS	3	20

Submental

One hole, flat backed for areas closer to the skin



Product Code	Diameter (mm)	Length (cm)
GI-210MTS	2	10
GI-315MTS	3	15

Open harvest

Guarded fully open tip, designed specifically for tissue harvest.



Product Code	Diameter (mm)	Length (cm)
GI-318HARV/TS	3	18