Advanced Lifecare

03 6343 4825

sales@advancedlifecare.com.au



RoMedic AutoSlide

Product Code: ALC218100

Getting into cars can be difficult and awkward, but once in, they are the key to freedom for many people. We have developed a device that makes it easier to get into cars and into a comfortable position. AutoSlide reduces the effort when turning and moving sideways to readjust the sitting position before and during the journey. If a little more assistance is needed, the handles at the sides of the seat cushion allow for giving a helping hand. The seat cushions are firmly attached to the backs of car seats, so users cannot slide forward. The device is very simple to fit and use - it's just a matter of starting the engine and driving off. The outer surface and straps are made of polyester. The...

Read more on our website

Features

- AutoSlide facilitates getting into and out of cars and reduces the effort when readjusting the sitting position before and during the journey.
- The soft, padded surface provides for a comfortable journey.
- The ULF sliding material reduces friction under buttocks.
- Robust handles enable carers to provide assistance.
- Secure fitting in the car seat prevents users from sliding forward.

355 Wellington Street, South Launceston.

Corner of Cridge St and Wellington St. Just near the pedestrian overpass.

advancedlifecare.com.au