

Enable Lifecare PTY LTD

ABN 86 620 684 586 ACN 620 684 586 support@enablelifecare.com.au
1300 370 370

enablelifecare.com.au

CONFIGURA COMFORT ASSESSMENT FORM

Prescriber's Name:			Client's Name:				
Number:			Add	lress:			
Email:			Nun	nher·			
Funding Body	/:						
Siz	Ze Areas highlighted in	n grey are standard specifica	ation				
Chair Size	Part Number	Seat Depth		Seat Height	Sea	t Width	Weight Capacity
		406 mm (16")	#1	388 mm (15") #4			
☐ SMALL	☐ CR5404 - Durate	ek 457 mm (18")		406 mm (16")	457	mm (18")	130 kg
	CR5422 - Vinyl	508 mm (20"))	431 mm (17") #2 457 mm (18") #3			
		406 mm (16")		439 mm (17") #4			
MEDIUM	CR5406 - Durate	ek		457 mm (18")	508 mm (20")		160 kg
	CR5424 - Vinyl	457 mm (18")		482 mm (19") #2			
		508 mm (20"))	508 mm (20") #3			
	CR5408 - Durate	406 mm (16")	#1	439 mm (17") #4			
LARGE		457 mm (18")	457 mm (18") 459 mm (18") 559 m		mm (22")	160 kg	
CR5426 - Vinyl		☐ 508 mm (20"))	508 mm (20") #3			
#1 Depth Adjustment #2 Hei		#2 Height Adjustment	<u></u>	^{#3} Height Adjustmer	nt 🗌	#4 Heigh	t Reduction 🗌
Kit 2" (50 m	Kit 1" (25 mm)	,	Kit 2" (50 mm)	44	Kit 1" (18	Tile AL	
CA2416		CA2421		CA2422			
5/12 / 12 5/12 /					9		
Pressure Posture							
Visco a Cushion (std) ☐ CushionAir+Pump CA2400 ☐				Waterfall Backrest (std) Laterals CA2413			
	0	999999	3				
			-			Со	me as each
					. –		







Enable Lifecare PTY LTD

ABN 86 620 684 586 ACN 620 684 586
support@enablelifecare.com.au
1300 370 370
enablelifecare.com.au

Transfer Options

Description	Long Description	Part Number	Select
Drop-down Armside Left-hand side (as sitting in the chair)	Transfers; Enables users to do a slideboard transfer in and out of the chair with ease. Can be bolted on at any time if the users needs change.	Small - Duratek CA2484 Small - Vinyl CA2492 Medium / Large - Duratek CA2488 Medium / Large - Vinyl CA2496	S
Drop-down Armside Right-hand side (as sitting in the chair)	Transfers; Enables users to do a slideboard transfer in and out of the chair with ease. Can be bolted on at any time if the users needs change.	Small - Duratek CA2486 Small - Vinyl CA2494 Medium / Large - Duratek CA2490 Medium / Large - Vinyl CA2498	S



Extras



Cognitive



Support

Armrests		2-1
Overhang on outside (std)		
Overhang on inside - reduces width between the armrests by 1.5" (38mm)		SC

2-Button Handset					
SC2316	Mary Common	Company			

Legrest Chan		
Small CA2520		
Medium CA2521		4
Large CA2522		9

Rise Setting		Legrest Setting	Handset Control	
Rise and Tilt (std) 23°	160kg SWL*	Standard	RHS (std)	
Medium Rise 14°	135kg SWL*	Standard + 1" (25 mm)	LHS	
Flat Rise 2°	110kg SWL*	Standard + 2" (50 mm)	Carer position	

^{*}MEDIUM & LARGE CHAIRS

Chair Assessment and Delivery Environmental Questionnaire

When carrying out an assessment please fill out the below form prior to ordering the chair.

Please give as much detail as possible.

ACCESS

ACCESS	
Is there parking/access for a large van?	
Are there steps to the property?	
Is it easy access to the property?	
INTERNAL	
Are the doors and corridors wide enough for a chair to go through? M	in. width for most chairs is 650mm.
Does the chair have to go upstairs? If yes, please specify number of s	steps and if a stair lift is fitted.
Is space being cleared for the chair?	
Is there a suitable power supply for the chair?	
INFECTION CONTROL	
Are there pets?	Is the house smoke free?

OTHER

Does the driver need to be aware of any other points? __