

RESPIRONICS

Sleep therapy

PAP retail accessories

Rediscover dreams

PAP accessories

to meet your patients' needs

Discover new ways to enhance your patients' therapy experience by offering a wide variety of convenient PAP accessories. Carrying these accessories can enable you to:

- Keep in contact with patients and build lasting relationships for future opportunities
- Help patients enjoy a better quality of life with a variety of accessories that allow them to conveniently travel with their sleep therapy system or store components discreetly at home
- Create a retail revenue opportunity and attract additional patients to your business





PAP battery kit

Give your patients the confidence and convenience of powering their PAP therapy device with the Philips Respironics battery kit. It's designed to provide an uninterruptible power supply (UPS) to the DreamStation and System One sleep therapy devices.

Designed for peace of mind

Your patients will know that when the power goes out, the battery pack's integrated uninterruptible power supply (UPS) will help reduce therapy disruptions. Continuously connecting the battery kit to the patient's Philips Respironics sleep therapy device's power supply keeps it charged and ready for travel or UPS protection.

Convenient for travel

The battery pack is our lightest and smallest ever, offering easy portability so patients can continue therapy away from home. It's great for camping, hotel use, or active duty military. The battery pack uses the power supply from the patient's Philips Respironics sleep therapy device to reduce the number of cords needed for travel. And it meets FAA requirements for airplane travel.

Benefits for patients

- The battery pack provides more than 14 hours of battery life* for PAP devices.
- It's compact and lightweight at just 2.25 lbs
- It meets FAA requirements (90W lithium ion battery pack
- No prescription is required.
- \cdot It comes with a six-month warrar
- Its lithium ion battery is environmentally conscious





*PAP run-time test conditions: CPAP mode, pressure 10 cm H2O, no heated humidification or heated tubing.

**PAP Lithium Ion Battery Kit for use with the Philips Respironics DreamStation and System One CPAPs, BiPAPs, autoSV, BiPAP AVAPS, and BiPAP S/T sleep therapy devices. Includes a rechargeable lithium ion battery pack, a PAP device cable, and a battery case. Carrying case

Components

The battery kit accommodates the Philips Respironics DreamStation and System One sleep therapy devices. It includes a rechargeable lithium ion battery pack, a PAP device cable, 2 adapters (for use with System One 50 series only), and a battery case.

Device cable

and adaptors

Supports DreamStation, System One 60 series, System One 50 series CPAPs and BiPAPs. Utilizes DreamStation or System One power supply. Additional power supply not included.

Specifications

Dimensions: 1.25" (3.175 cm) x 5.25" (13.335 cm) x 9.0" (22.86 cm) Weight: 2.25 lbs. / 1.02 kg

Ordering information	Part number
PAP Lithium Ion Battery Kit**	1124009
Replacement Parts	Part number
Li Ion battery, PAP device cable	1121736
Li Ion Battery, Cable Adapter (50 series adapter, qty 2)	1127799
Li Ion Battery Case	1126947
DreamStation 65W power supply	1120136
DreamStation 80W power supply	1118499
DreamStation/System One power cord	1005894
System One 60W power supply (for H and HS models)	1091398

System One 80W power supply (for TS model only) 1091399



DreamStation travel case

The DreamStation travel case* gives patients a convenient, compact, and discreet way to travel with their DreamStation sleep therapy device and supplies.

Convenient

The case includes a DreamStation 65W power supply that is lighter and smaller than the DreamStation 80W power supply. The 65W power supply powers the DreamStation** device, and the included USB port can be used to charge items such as a smartphone and other electronic devices during travel.

Compact

The lightweight, hard-fabric case offers an ideal combination of protection and portability. Its individual compartments accommodate the DreamStation** device and neatly organize the supplies. An included mesh pouch holds the tubing securely in place.

Discreet

The case is designed to make traveling with sleep therapy equipment less cumbersome. A molded insert liner helps the patient thoroughly pack the PAP device and accessories. The case can be attached to a roller suitcase and its zip top facilitates airport security screening. On a plane, it easily fits in the overhead bin or under the seat.

Benefits for patients

- mind that the sleep therapy device will be

Luggage strap

Detachable shoulder strap

Power supply with **USB** port DreamStation Tubing

Power supply

cord

Components

without humidifier

Mask

The case accommodates the Philips Respironics DreamStation sleep therapy device (blower only, no humidifier), the patient's mask and tubing, a 65W power supply w/ USB port, and a removable luggage strap.

Specifications

Dimensions: 10.42" (26.46 cm) x 3.94" (10 cm) x 13.78" (35 cm)

Ordering information	Part number
Travel case, 65W power supply w/ USB port	1120135
Travel case, 65W power supply w/ USB port, 5/case	1124301

*DreamStation, mask, and tubing not included

**Without humidifier



PAP travel briefcase with detachable laptop bag

The PAP travel briefcase* gives patients an all-inone carrying option that accommodates most PAP therapy systems and humidifiers, a laptop, and other small items. It makes it easy and efficient to travel with a sleep therapy system.

The travel briefcase features two bags that conveniently zip together into one carry-on. The therapy device and accessories fit snugly into customizable compartments on one side of the case, while the laptop and its accessories fit securely into the ample-sized pockets of the other side. There's even space for books, magazines, and other small travel items.

Benefits for patients

- The case's compartments protect the sleep therapy device and accessories during travel.
- Its zip top facilitates airport security screening
- It fits within the FAA maximum carry-on size
- The lightweight case offers easy portability
- It's designed to attach to a roller suitcase.



Detachable shoulder strap

Removable laptop bag **Power supply**

Tubing

DreamStation with humidifier

Components

The briefcase accommodates most PAP therapy devices and accessories and standard-size laptop computers.

Specifications

Dimensions: 16.14" (40.99 cm) x 13.58" (34.49 cm) x 7.87" (19.99 cm) Color: Black

Ordering information	Part number
PAP travel briefcase	1114784

*DreamStation, power supply, and tubing not included



Mask and tubing bedside organizer

The bedside organizer gives patients a way to conveniently store their mask and tubing without cluttering a nightstand – both at home and during travel. It also keeps the components dust-free and ready for use.

The organizer has a weighted flap to allow it to hang over a nightstand or between a patient's mattress and box spring. It even features a zippered compartment for the mask and tubing, making the organizer a great travel companion.

Benefits for patients

- The organizer provides a clean, discreet place to store the mask and tubing while not in use.
- Its design looks right at home in the bedroom
- · It can be used on a nightstand or in between a mattress and box spring.

Specifications

Dimensions: 7.99" (20.3 cm) x 10.94" (27.8 cm) x 3.74" (9.5 cm) Color: Light gray

Ordering information	Part number
Bedside organizer	1122135
Bedside organizer 10/case	1124224





LiquiCell[®] nasal CPAP cushions

LiquiCell* is a thin, pliable cushion infused with fluid that can be applied virtually anywhere on the face, making it versatile. The cushion enhances comfort by acting as a barrier between the user's face and the mask. Gentle adhesive keeps the cushion securely in place.

Part number

1109097

Ordering information

LiquiCell nasal CPAP cushions, pack of 30

*LiquiCell is a registered trademark of Birchwood Laboratories, Inc.

Order your PAP accessories today

For more information, call 1-800-345-6443 or visit www.philips.com/sleeptherapy.



© 2016 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice.

www.philips.com/respironics

Caution: US federal law restricts these devices to sale by or on the order of a physician.

RRDPGH EB 7/25/16 MCI 4107650 PN 1129986 1010 Murry Ridge Lane, Murrysville, PA 15668 800-345-6443 • 724-387-4000