



SparMED Aps  
 Ryttermarken 2  
 3520 Farum  
 Denmark



140 Hale Street  
 Haverhill, MA 01830  
[qc@embryotech.com](mailto:qc@embryotech.com)

**ELI Accession Number:** S2806-0117SPAR

**Date of completion:** 01-06-2017

**Lot numbers:** 60465, 80350, 70355

**Reference numbers:** OODIH, OODH, OODSF

**Description of test articles:** Oosafe CO<sub>2</sub> Disinfection, Oosafe Hand Disinfection, Oosafe Surface Disinfection

**Assay system requested by customer:** An incubator was cleaned with the test article (Oosafe CO<sub>2</sub> Disinfection). Post cleaning 1mL of each test articles (Oosafe Hand Disinfection and Oosafe Surface Disinfection) was placed in a Petri dish and placed in the cleaned incubator. A prepared sperm sample was then placed in the cleaned incubator with the test articles for 24-hours. The forward progressive motility was read and recorded at 24-hours.

**Results:**

	<u>Initial motility</u>	<u>24hour motility</u>
<b>Test Article</b> (incubator ELI-182)	100 %	99 %
<b>Control</b> (incubator ELI-243)	100 %	100 %

**Summary of observations:** The motility of the sperm remained consistent in the incubator ELI-182 (that was cleaned with and that contained the test articles) and the incubator ELI-243 (the control sperm).

  
 \_\_\_\_\_  
 Signature  
 Study Director

01-06-2017  
 \_\_\_\_\_  
 Date

  
 \_\_\_\_\_  
 Signature  
 Quality Reviewer

01-06-2017  
 \_\_\_\_\_  
 Date