



SparMED Aps  
Ryttermarken 2  
3520 Farum  
Denmark



140 Hale Street  
Haverhill, MA 01830  
[qc@embryotech.com](mailto:qc@embryotech.com)

**ELI Accession Number:** S2712-1016SPAR    **Date of completion:** 10-18-2016

**Lot number:** 63504

**Reference number:** OODW-70

**Description of test article:** Oosafe® Disinfection Wipes

**Assay system requested by customer:** An incubator was cleaned with the test article. Post cleaning a prepared sperm sample was placed in the cleaned incubator for 24-hours. The forward progressive motility was read and recorded at 24-hours.

**Results:**

	<u>Initial motility</u>	<u>24hour motility</u>
<b>Test Article</b> (incubator ELI-038)	96 %	94 %
<b>Control</b> (incubator ELI-068)	96 %	96 %

**Summary of observations:** All test and control sperm was prepared from the same donor. The incubator ELI-068 (the control sperm) had a 96% forward progressive motility at 24-hours. The incubator ELI-038 (that was cleaned with the test article) had a 94% forward progressive motility at 24-hours.

Signature  
Study Director

10-19-2016  
Date

Signature  
Quality Reviewer

10-19-2016  
Date