

# ICONOCLAST®

## EQUINE SUPPORT BOOTS

### MEASUREMENT & SIZING CHART

The sizing chart below is appropriate for all Iconoclast® Orthopedic and Rehabilitation Boots.

We recommend using a soft tape measure. If you do not have a soft tape measure you can use something like string, wrap it around the fetlock, mark where the end meets the rest of the string and then simply lay it flat next to a tape measure, ruler or yard stick to get your measurement.

Be sure to measure both the FRONT & BACK fetlocks as they are often different sizes. Customers regularly order boots in one size for the front legs and a different size for the hind legs.



Small: 8.5 - 9.5 in. / 21.59 - 24.13 centimeters  
Medium: 9.5 - 10.5 in. / 24.13 - 26.67 centimeters  
Large: 10.5 - 11.5 in. / 26.67 - 29.21 centimeters  
XL: 11.5 - 12.5 in. / 29.21 - 31.75 centimeters  
XXL: 12.5 - 13.5 in. / 31.75 - 34.29 centimeters