



Medical Device Details

Sponsor: CARSL Consulting
P O Box 766
Hastings
New Zealand

WAND Reference: 180615-WAND-6QIMAA

Sponsors Own Reference: Precise Medical Supplies - Cannula, liposuction


GMDN: Liposuction system cannula, single-use [36107]

Class: Ila

Intended Purpose: A rigid tube that is inserted percutaneously under the skin, and is used during the removal of fat deposits. An appropriate unit is used to provide the suction. This cannula is single use.

Manufacturer: INEX Z.A la Gobette
Puisseux le Hauberger
FRANCE

ARTG ID: 157914

 > The device is supplied sterile.

- > Sterilisation method: Ethylene oxide
- > The device is intended to be invasive.
- > The device is intended for single use.
- > The device is not an active device.
- > The device does not contain material or ingredients of microbial origin.
- > The device does not contain material or ingredients of recombinant origin.
- > The device does not contain material or ingredients manufactured or formulated using a genetically modified organism.
- > The device does not contain material or ingredients of human origin.
- > The device does not contain human blood or its components.
- > The device consists of: Single product
- > The device does not contain material or ingredients of animal origin.
- > The device is not medicated.
- > The device is not formulated.
- > The product does not contain a medicine that has consent for marketing in New Zealand.
- > The product does not contain a medical device which incorporates a medicine as an integral part and that has an action ancillary to the device.

Status: Active