



## Medical Device Details

Sponsor: CARSL Consulting  
P O Box 766  
Hastings  
New Zealand

WAND Reference: 180615-WAND-6QIMJJ

Sponsors Own Reference: Precise Medical Supplies - Needle, hypodermic


GMDN: Hypodermic needle, single-use [12745]

Class: Ila

Intended Purpose: To administer fluids to a patient or to withdraw fluids from a patient

Manufacturer: TSK Laboratory Japan 2-1-5 Hirayanagi-cho  
Tochigi-Shi  
Tochigi-Ken  
JAPAN

ARTG ID: 194674

-  > The device is supplied sterile.
- > Sterilisation method: Gamma irradiation
  - > The device is intended to be invasive.
  - > The device is intended for single use.
  - > The device is not an active device.
  - > The device does not contain material or ingredients of microbial origin.
  - > The device does not contain material or ingredients of recombinant origin.
  - > The device does not contain material or ingredients manufactured or formulated using a genetically modified organism.
  - > The device does not contain material or ingredients of human origin.
  - > The device does not contain human blood or its components.
  - > The device consists of: Single product
  - > The device does not contain material or ingredients of animal origin.
  - > The device is not medicated.
  - > The device is not formulated.
  - > The product does not contain a medicine that has consent for marketing in New Zealand.
  - > The product does not contain a medical device which incorporates a medicine as an integral part and that has an action ancillary to the device.

Status: Active