

Zircon Plus



Product Features:

Zircon Plus Sintering Furnace delivers delicacy and beauty to your zirconia. It features molybdenum disilicide heater that provides clean, reliable, and consistent results. It is also user-friendly with very simple operation. The unit is applicable to various types of zirconia materials. Large furnace allows full mouth case

Product Specifications:

Power Requirement:

AC 220V. 50/60Hz. Max. 2.5kVA

Dimension:

17.8 (w) x 28.6 (h) x 21.3 (d) inches

Display:

3.5inch LCD Color Display

Firing Platform:

Ø95mm (3.74 inch)

Weight:

110lbs (50kg)

Max Temperature:

1,600°C (2,912°F)

Programs:

100 (20 Steps)

SD Card Slot:

Software Update

Program Back Up