


Walking sling

New and improved sling for active users who need extra support and safety during rehabilitation exercises



Usability

- For users who need support with walking or standing training.
- For training with Ergo Trainer or stationary hoists.
- Strap system on the sling can be applied to most 2-point spreader bars.



Features

- The sling is allowing many customising options. It is possible to adjust the spread around the abdomen, the length of the leg straps and shoulder straps.
- The sling ensures safety and support, in case the user lost balance under training.



Walking sling has

- removable padding on the shoulder straps.
- removable padding on the leg straps.
- elastic back, ensuring a good fit for the user.
- coloured strap on upper back with size and mid-line mark, which ensure correct placement of the sling.
- **NEW!** Elastic sides, ensuring a better fit for the user - almost with the same function as a kidney belt.
- **NEW!** The leg strap buckles are sewn on in an oblique angle to make the strap fit better.
- **NEW!** The design of the sling has been changed, making it possible to place and buckle it from behind. It is also still possible to buckle it on the front.



Item no.	Size	Model
009-09572	XXS	Walking sling
009-09883	XS	Walking sling
009-09884	S	Walking sling
009-09885	M	Walking sling
009-09886	L	Walking sling
009-09887	XL	Walking sling



Sling size guide (cm)

Colour code is presented on the size chart.
The image below shows red colour code, which means that the sling is in size S (Small).



Technical information

SWL 180kg



	XXS	XS	S	M	L	XL
Belly dimension	62 – 75 cm	76 – 84 cm	85 – 92 cm	93 – 105 cm	106 - 114 cm	116 - 130 cm

See the user manual for detailed specifications. All figures and measurements are for guidance only and may differ in actual situations. Ergolet reserves the right to changes without notice.

Product photos may include accessories. Consult your dealer.

WinnCare Nordic ApS
Taarnborgvej 12C
4220 Korsoer
Denmark
tel. +45 70 27 37 20
e-mail: info@winnCare.dk

WINNCARE
N o r d i c