

Enertor® 3D Print Bespoke Prescription Form

Tel. [01473] 834 891 | www.medical.enertor.com | E-mail: Orders@enertor.com

I agree with Enertor® Terms and Conditions

Standard [12 days]



ENERTOR
MEDICAL
Admin Use Only

Practitioner/Company

Address Post Code

Phone Email

Patient Name

No. Of Pairs Repeat pair Shoe Size Weight Kg

Attachments: Template Foam Box Cast 3D Scan Footscan Footwear Other

Order #

Received

Despatched

Account

Clinical Support

1. Shell Details

Shell Style:

- Classic*
- Heel Pain
- Slim Fit
- Fashion Fit
- Marathon
- Cobra
- Football
- Rugby

Shell Thickness:

- [2mm] Flexible
- [3mm] Semi-Flexible*
- [4mm] Rigid

Heel Cup Depth:

- [10mm] Shallow
- [15mm] Deep*
- [20mm] Extra Deep

2. Top Cover

Length: Shell Length Sulcus Length Full Length* cm

Width: Standard Fit* Narrow Fit Wide Fit cm

Top Layer:

- ELITE*
- VILEDON¹
- SMART LEATHER¹
- SMART SUEDE¹

Cushion Layer²:

- PX1® 3mm*
- EXTRAFOAM SOFT 3mm
- PORON 4000: 1.6mm 3.2mm 6.4mm
- PORON94: 1.6mm 3.2mm 6.4mm

EVA³:

- AQUA 2mm
- AQUA 3mm
- LUNASOFT 3mm
- HARD EVA 3mm

*USED AS STANDARD IF NOTHING CHOSEN, ¹ADDITIONAL CHARGE APPLIES,
²MUST BE USED WITH TOP LAYER, ³CAN BE USED WITH OR WITHOUT TOP LAYER.

3. Profile and Posting

Shell Arch Height Left Right

High [+3mm]

Normal* [To Cast]

Low [-3mm]

Rearfoot Wedge

Left Right

Medial°

Lateral°

No Heel* 3D Printed Heel EVA Heel

Forefoot Wedge

Left Right

Medial°

Lateral°

Built In Shell* 3D Printed EVA

4. Shell Modification

	Left	Right
Medial Heel Flange	<input type="checkbox"/>	<input type="checkbox"/>
Lateral Heel Flange	<input type="checkbox"/>	<input type="checkbox"/>
Medial Midfoot Flange	<input type="checkbox"/>	<input type="checkbox"/>
Lateral Midfoot Flange	<input type="checkbox"/>	<input type="checkbox"/>
Medial Heel Skive [mm]	<input type="text"/>	<input type="text"/>
Lateral Heel Skive [mm]	<input type="text"/>	<input type="text"/>
Heel Raise [mm]	<input type="text"/>	<input type="text"/>
1st Ray Cut Out	<input type="checkbox"/>	<input type="checkbox"/>
5th Ray Cut Out	<input type="checkbox"/>	<input type="checkbox"/>
1st MTPJ Cut Out	<input type="checkbox"/>	<input type="checkbox"/>
5th MTPJ Cut Out	<input type="checkbox"/>	<input type="checkbox"/>
PX1® Underfill	<input type="checkbox"/>	<input type="checkbox"/>
PORON Underfill	<input type="checkbox"/>	<input type="checkbox"/>

5. Padding

	Left [mm]	Right [mm]	PORON*	EVA ⁴	ES ⁵	PX1® ⁶
Medial Arch Pad	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Heel Cushion Pad	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Heel Spur Pad	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Morton's Extension	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Reverse Mortons Extension	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Neuroma Pad	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	5 4 3 2 1	1 2 3 4 5				
Neuroma Pad Location	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>				
Plantar Metatarsal Pad [PMP]	<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	5 4 3 2 1	1 2 3 4 5				
PMP Location	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	- LEAVE BLANK FOR CUTOUTS			

⁴WHITE EVA AND ⁵EXTRAFOAM SOFT - AVAILABLE IN VARIOUS THICKNESSES.

⁶PX1® AVAILABLE FROM 3mm, *ALL PADDING IN 1.6mm PORON UNLESS STATED OTHERWISE.

Additional Instructions and drawings:

