



Enertor® 3D Print Prescription Form

Tel. 01473 834891 www.medical.enertor.com Email: orders@enertor.com

I agree with Enertor® Terms and Conditions

12 working day turnaround as standard.

Admin Use Only	
Order #	
Received Date	
Despatch Date	
Account #	
<i>Clinical support completed</i>	

Practitioner							Post Code		
Address									
Phone									
Patient Name									
No. of Pairs	Repeat Pair	Enclosed -	3D Scan	Foam Box	POP Cast	Footscan Files	Template	Footwear	
			Email	Shoe Size	Weight	Kg			

Heel Cup Option		Shell Thickness		Top Cover Options	
Shallow Heel Cup (1cm)		2mm (flexible)		PX1® (with Topflight Fabric) */**	
(Deep Heel Cup** (1.5cm)		3mm (semi flexible)		PX1® & Elite	
Extra Deep Heel Cup (2cm)		4mm (rigid)		Ultra Cushioning (Beige EVA)	
				Cushioning (Grey EVA)	
				Cushioning (Aqua EVA)	
				Support (Black EVA)	
				Lunasoft (water resistant EVA)	
				Elite & Poron (1.6mm poron)	
				Elite (1mm Faux Leather)	
				Custom Elite (If purchased)	
				Custom Elite & PX1® (If purchased)	
				Custom Elite & Poron (If purchased)	

3D Print Style	3D Print Sport Styles
Classic**	Marathon
Heel Pain	Football
Slimfit	Golf
Fashion Fit	
Cobra	

*Please Note: ** text denotes orthotic style and materials used on if none selected.*

Top Cover Length	¼ Length	Sulcus Length	Full Length*
Shoe width	Standard Fit*	Narrow Fit	Wide Fit

2. Device Profile and Posting			Rearfoot Posting			Forefoot Posting		
Arch Height	Left	Right	Intrinsic	Extrinsic	Intrinsic	Extrinsic	Intrinsic	Extrinsic
Low (4mm lower than cast)			Medial Left deg.	Medial Right deg.	Medial Left deg.	Medial Right deg.		
Normal (2mm lower than cast)			Lateral Left deg.	Lateral Right deg.	Lateral Left deg.	Lateral Right deg.		
High (to cast)			Post Material: All heels will be 3D printed			Post Material 3D printed EVA		

3. Shell Modifications		4. Padding	
Left	Right	Left (mm)	Right (mm) PX1®
Medial Heel Flange		Medial Arch Pad	
Lateral Heel Flange		Heel Cushion Pad	
Medial Midfoot Flange		Heel Spur Pad	
Lateral Midfoot Flange		Morton's Extension EVA Poron	
Medial Heel Skive	mm mm	Reverse Morton's Ext EVA Poron	
Lateral Heel Skive	mm mm	Neuroma Pad	
Heel Raise	mm mm		
1st Ray cut out			5 4 3 2 1 1 2 3 4 5
5th Ray cut out		Neuroma Pad Location	
1st MTPJ cut out		Plantar Metatarsal Pad(PMP)	
5th MTPJ cut out		PMP Cut Outs	5 4 3 2 1 1 2 3 4 5
Underfill (Carbon only) (PX1® used as standard)			

NB: PX1® available from 3mm, all padding in 1.6mm Poron unless stated otherwise.

Additional Instructions:

