



# Enertor® Bespoke Prescription Form

Tel. 01473 834891 [www.rightstride.com](http://www.rightstride.com) Email: [orders@enertor.com](mailto:orders@enertor.com)

## Admin Use Only

Order #  
Received Date  
Despatch Date  
Account #

I agree with Enertor® Terms and Conditions

Standard (8 days)    Speedy (3 Days)    **Express (24Hrs)**

Practitioner

Address

Phone

Patient Name

No. of Pairs

Repeat Pair

Enclosed -

3D Scan

Foam Box

POP Cast

Footscan Files

Template

Footwear

Post Code

Clinical support completed

Email

Shoe Size

Weight

Kg

## 1. Orthotic Style and Materials

### EVA Orthotic Style

Classic\*

Sport/Activity

EVA Base Option    Full    3/4

Ultra Soft (Beige)

Soft (Blue)

Supportive (White)\*

Maximum Support (Black)

Please note - This option determines the length of the milled EVA section of the orthotic.

Please Note: \* text denotes orthotic style and materials used on EVA if none selected.

### Carbon Base Option

Shallow Heel Cup (1cm)

(Deep Heel Cup\*\* (1.5cm)

Extra Deep Heel Cup (2cm)

### Carbon Resin Orthotic

Classic\*\*

Heel Pain

Slimfit

Fashion Fit

Cobra

### Carbon Sport Styles

Marathon

Football

Golf

### Shell Thickness

2mm (Flexible)

3mm (Semi-flexible)

4mm (Rigid)

Please Note: \*\* text denotes orthotic style and materials used on if none selected.

### Top Cover Options

PX1® (with Topflight Fabric) \*/\*\*

PX1® & Elite

Ultra Cushioning (Beige EVA)

Cushioning (Grey EVA)

Cushioning (Aqua EVA)

Support (Black EVA)

Lunasoft (water resistant EVA)

Elite & Poron (1.6mm poron)

Elite (1mm Faux Leather)

Custom Elite (If purchased)

Custom Elite & PX1® (If purchased)

Custom Elite & Poron (If purchased)

Top Cover Length    ¾ Length

Sulcus Length

Full Length\*

Shoe width    Standard Fit\*

Narrow Fit

Wide Fit

## 2. Device Profile and Posting

Arch Height

Left

Right

LOW (4mm lower than cast)

Normal (2mm lower than

cast) High (to cast)

### Rearfoot Posting

Intrinsic

Extrinsic

Medial Left

deg.

Medial Right

deg.

Lateral Left

deg.

Lateral Right

deg.

### Stabiliser Material

Direct Milled

EVA

### Forefoot Posting

Intrinsic

Extrinsic

Medial Left

deg.

Medial Right

deg.

Lateral Left

deg.

Lateral Right

deg.

## 3. Shell Modifications

Medial Heel Flange (Carbon only)

Lateral Heel Flange (Carbon only)

Medial Midfoot Flange (Carbon only)

Lateral Midfoot Flange (Carbon only)

Medial Heel Skive

mm

mm

Lateral Heel Skive

mm

mm

Heel Raise

mm

mm

1st Ray cut out

5th Ray cut out

1st MTPJ cut out

5th MTPJ cut out

Underfill (Carbon only) (PX1® used as standard)

## 4. Padding

Medial Arch Pad

Heel Cushion Pad

Heel Spur Pad

Morton's Extension

EVA

Poron

Reverse Morton's Ext

EVA

Poron

Neuroma Pad

5 4 3 2 1

1 2 3 4 5

Neuroma Pad Location

Plantar Metatarsal Pad(PMP)

5 4 3 2 1

1 2 3 4 5

PMP Cut Outs

NB: PX1® available from 3mm, all padding in 1.6mm Poron unless stated otherwise.

## Additional Instructions:

