



Economy Automatic Blood Pressure Monitor, Upper Arm

Item #	Product Description	UOM
BP3600	Economy Automatic Blood Pressure Monitor, Upper Arm	1/ea, 24/carton

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory Stores 120 Readings
- Average Last 3 Results

specifications	
Unit Dimensions	5.51" x 3.86" x 1.89"
Carton Size	18.9" x 17.7" x 17.7"
Measurement Method	Oscillometric Method
Measurement Range	Pressure: 0mmHg~300 mmHg Pressure: +/- 3 mmHg
Power Supply	4x AAA Batteries or AC Adapter
Auto Power Off	3 Minutes After Last Reading
Cuff Size	8.5"-14"
Memory Capacity	120 Readings with Date & Time
Operation Environment	50°F to 104°F, <85% RH
LCD Size	2.17" x 3.31"
Warranty	2 Years

accessories		
BP2000	Plus Size Blood Pressure Cuff For BP 2200 and BP2400	1/ea
BP2001	Large, 8.5" x 14"	1/ea
BP2610	A/C Adapter for BP2600 and BP2400	1/ea