

Pro-Vent™ Mask Configurations

Configuration	Part #	NDC #	HCPCS #
Bubbles the Fish™ II (Ages 3+ years)	#044F7248	#83490-0440-68	A7015
PARI LC® Adult Mask	#044F7251	#83490-0440-78	A7015


Features

- Front-load design directs aerosol toward the patient's mouth
- Comfortable aerosol delivery with soft, flexible masks
- Pediatric character design makes treatments fun for kids!



The Difference is in the Delivery

■ Bubbles™ II with PARI LC® Plus
 ■ Salter Mask with Misty Neb



Accessories:

PARI LC® Sprint = **023F35**
 PARI LC® Plus = **022F81**
 PARI Pals™ Activity Book = **044D0029**
 PARI Pals™ Stickers = **044D0030**

KEY FINDINGS: A higher Inhaled Mass combined with lower facial and eye deposition gives patients increased drug delivery to the lungs with reduced exposure to the eyes and skin.

¹Facial and Ocular Deposition of Nebulized Budesonide. K. Harris, G. Smaldone. CHEST 2008; 133:482-488. Tables 2-4.
²All registered marks and trademarks belong to their respective companies.

PRO Vent™
 TECHNOLOGY



Fun, safe, optimized delivery of aerosol.

Only the PARI adult mask and Bubbles The Fish™ II pediatric mask are designed with Pro-Vent™ technology, proven to deliver more medication while reducing facial and eye deposition.

From the leader in clinically proven aerosol delivery