

Soft Nasal Cannula - Non-Kink Tubing



Made of extra-soft disposable material for long-term use. Comes standard with non-kink tubing in various lengths up to 50 feet. Adjustable tubing fits comfortably over the ears to hold the cannula firmly in place, affords complete freedom of movement.

Adult available in all lengths; Pediatric available in 7-foot length.

Curved

| | |
|------------|---------------------|
| SOFT 001 | Cannula Only, 50/cs |
| SOFT 004 | 4' Length, 50/cs |
| SOFT 007 | 7' Length, 50/cs |
| SOFT 007 P | 7' Length, 50/cs |
| SOFT 014 | 14' Length, 50/cs |
| SOFT 025 | 25' Length, 25/cs |
| SOFT 050 | 50' Length, 20/cs |

COMPARE TO:

| | |
|---|---|
|  |  |
| SOFT 007 | 1110 |
| SOFT 014 | 1810 |
| SOFT 025 | 1812 |
| SOFT 050 | 1814 |



Curved SOFT 007



Straight SOFT 107



Compares to Salter™
Cozy Cannula SOFT 207


Straight Tipped

| | |
|------------|---------------------|
| SOFT 100 | Cannula Only, 50/cs |
| SOFT 104 | 4' Length, 50/cs |
| SOFT 107 | 7' Length, 50/cs |
| SOFT 107 P | 7' Length, 50/cs |
| SOFT 114 | 14' Length, 50/cs |
| SOFT 125 | 25' Length, 25/cs |
| SOFT 150 | 50' Length, 20/cs |

The Cozy Cannula, Non Kinking Tubing

| | |
|------------|-------------------|
| SOFT 200 | Cannula, 50/cs |
| SOFT 204 | 4' Length, 50/cs |
| SOFT 207 | 7' Length, 50/cs |
| SOFT 207 P | 7' Length, 50/cs |
| SOFT 214 | 14' Length, 50/cs |
| SOFT 225 | 25' Length, 25/cs |
| SOFT 250 | 50' Length, 20/cs |

COMPARE TO:

| | | |
|---|---|--|
|  |  |  |
| SOFT 200 | 1821 | 1606 |
| SOFT 204 | 1827 | 1600-4 |
| SOFT 207 | 1820 | 1600 |
| SOFT 214 | 1817 | 1600-14 |
| SOFT 225 | 1823 | 1600-25 |
| SOFT 250 | 1831 | 1600-50 |



Dual Lumen
Consaver Cannula SOFT307

Dual Lumen

Consaver Cannula

| | |
|----------|------------------|
| SOFT 304 | 4' Length, 25/cs |
| SOFT 307 | 7' Length, 25/cs |

Fits All Cannulas with Universal Connector

Fits standard barbed nipples and DISS (Diameter Indexed Stainless Steel) fittings. Comes in straight tipped style Adult & Pediatric cannula's only.

| | |
|------------|------------------------------|
| SOFT 107U | 7' Cannula, Adult, 50/cs |
| SOFT 107PU | 7' Cannula, Pediatric, 50/cs |



Universal Connector

All of our tubing is made in translucent green to be in compliance with JACO/OSHA tripping hazard recommendations.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*