

DeVilbiss Traveler® Portable Compressor Nebulizer System



DeVilbiss Traveler Portable Compressor Nebulizer Unit

The **perfect** travel companion.

The **Traveler Portable Compressor Nebulizer System** is designed to help active customers get more out of their aerosol therapy, and **more out of life**. It's never been easier or more convenient to nebulize outside of the home.

System includes compressor unit, disposable and reusable nebulizers, **five extra filters**, tubing adapter, adult and **pediatric** masks, AC adapter, DC adapter, carrying case, instruction guide and **instructional DVD**. Units can be ordered with **optional battery**.

Only DeVilbiss Traveler offers both adult and **pediatric masks**, a **disposable** and reusable nebulizer, an instructional **DVD** and a **full-year** battery warranty.

Performance Features

- Delivers a **consistent and concentrated treatment** of small particles that penetrates deeply into the lungs, delivering medication where it's needed most
- **Fast treatment time** (approximately 10 minutes) with DeVilbiss 800 Series Reusable Nebulizer encourages patient compliance
- Can be used with **any manufacturer's nebulizer** – treatment times will vary
- **45-minute battery life** far outlasts treatment time
- The 4-hour battery charge time is the **shortest in the category**

Design Features

- Easy, one-touch operation
- Compact and lightweight
- AC adapter and DC car adapter
- One-year battery and three-year compressor warranty

Model 6910D-DR Includes:

- Compressor unit
- Reusable nebulizer with tubing
- Disposable nebulizer
- 5 replacement air filters and filter door
- Tubing adapter
- Adult and pediatric masks
- AC adapter
- DC car adapter
- Easy tote carrying case
- Instructional manual and DVD



DeVilbiss Traveler Portable
Compressor Nebulizer System
Model 6910P-DR

This is also available with a battery as Model 6910P-DR and comes with all of the above as well as a rechargeable battery pack.

Free Consumer Marketing Program to Generate Sales

- **Comprehensive Traveler marketing kit** includes counter card with reminder tear pad, product flyer (stapled to prescription bags) and patient postcard mailers to **generate business**

Technical Specifications

MMAD	5 microns	Sound Level	54 dBA
Size w/o Battery	3.34" W x 2.2" H x 3.8" D	Max. Pressure	30 psig or greater
Size w/ Battery	3.34" W x 2.2" H x 5.3" D	Free Air Flow	6 lpm or greater
Weight w/o Battery	.8 lbs	Intermittent Use	30 minutes on/30 minutes off
Weight w/ Battery	1.5 lbs	Carton Size	8 units
Electrical	100-240V, 50-60Hz, 15VDC, 1.5A	Battery Life	45 minutes



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com