

NEW!

**PATENTED
MOTION
SENSING!**



OM-900M

**THE EVOLUTION™ MOTION OXYGEN CONSERVER
WILL SENSE PATIENT ACTIVITY AND TURN OXYGEN FLOW
UP AND DOWN AS NEEDED, WITHIN SECONDS**

Auto-adjusting conservers are the future of oxygen therapy, optimizing both patient saturation and maximum conservation. Satisfy your patients' and referrals' wishes while decreasing your delivery costs, all with one unit.



Evolution™ Motion, Electronic Oxygen Conserving

1P OM-900M 1/ea

FEATURES:

- Programmable Rest and Active setting equivalent to 1-6 liters per minute.
- Delivers 5:1 savings ratio at all settings.
- Minimum 1 year battery life utilizing 2 AA Batteries (avg. 4 hrs. per day usage).
- Delivers a uniform oxygen pulse from 14-40 breaths per minute.
- Liter flow increments: 1, 2, 3, 4, 5, 6 Plus CF preset @ 2 LPM.
- Operating pressure range: 200-3000 PSI.
- Senses a breath/delivers oxygen in the first half of the inspiratory cycle.
- Weighs only 14.9 ounces with batteries.
- Auto shutoff feature conserves battery life.
- Fits easily in virtually any carry bag.
- Limited 2 Year Warranty.

**THE EVOLUTION™ MOTION OXYGEN CONSERVER WILL SENSE PATIENT
ACTIVITY AND TURN OXYGEN UP AND DOWN AS NEEDED WITHIN SECONDS.**

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

1P = 1-Pack Available