

Briutcare[®]

Changing the way you breathe

New Generation Smart Mesh Nebulizer



Super-Silent Effect

Working noise less than 35 dB.
Automatically shutdown when no liquid detected.
< 40dB, even quieter than whisper



Rechargeable

Embedded lithium battery with USB charging.
4 hours' work



One-Button Operation

Simple and easy to use for all consumers.
Adjustable nebulization rate
0.3ml/min, 0.4ml/min, 0.5ml/min



Powerful & Portable

4-hour long nebulization.
Pocket size, easy to carry around.
100g, lighter than a smartphone



Patented I.A.D Technology

We've created a smart respiratory sensor which allows the nebulizer to atomize medicine while inhaling, and stop atomizing while exhaling.



App Game Guided Nebulization

Children love games but usually don't enjoy the nebulization process, so we created ours to assist parents and children improve the nebulization experience!



Ultrafine Particles

Thanks to our patented 'Vibration Mesh Technology', 80% of particles are smaller than 5µm, with a MMAD of only 2.2µm. Thus, less dosage time with better efficiency is realized!



World's first personalized nebulization rate available!

We make it possible for patients to adjust the nebulization rate and choose the best model and treatment for themselves.



*Atomizing Cup is recommended to be replaced every 2-3 months