




HCPCS
E0570



Compartment Style Compressor Nebulizer

Item #	Product Description	UOM
MQ6000	Compartment Style Compressor Nebulizer	1/ea

- Includes disposable neb kit and tubing
- Attractive, Sleek Design
- Durable, Quiet System
- Piston-Driven Compressor
- Convenient Storage Compartment

specifications	
Liter Flow	≥ 8 lpm
Nebulizer Medication Capacity	5ml
Particle Size	0.5 µm to 5.0 µm
Maximum Pressure	≥ 30 psi
Operating Pressure	≥ 10 psi
Operating Temperature Range	10° to 40°C (50° to 104°F)
Operating Humidity Range	10% to 95% RH
Storage Temperature Range	-20° to 70°C (-4° to 158° F)
Storage Humidity Range	10% to 95% RH
Power Requirements	120V, 50-60 Hz
Power Consumption	80 VA
Noise Level	55 dBA
Weight	4.3 lb
Warranty	Limited 5 Year

accessories		UOM
NEB KIT 500	Disposable Neb Kit with Mouthpiece, "T," and 7' Tubing	50/cs
MASK 001A	Adult Aerosol Mask	50/cs
18254	Reusable Neb Kit with Mouthpiece, "T," and 7' Tubing	10/cs