

## Oxygen Tubing - Non Kink


Made of soft disposable material. Three channel safety construction provides uninterrupted flow of oxygen. When uncoiled, tubing lies out straight...it has no memory. Thus, "torque" or pull on the cannula is eliminated. Therefore, tubing will not coil up on the floor.

TUB NK 07	7' Length, 50/cs
TUB NK 14	14' Length, 50/cs
TUB NK 25	25' Length, 25/cs
TUB NK 40	40' Length, 20/cs
TUB NK 50	50' Length, 20/cs

COMPARE TO:		
	MEDX <small>A Division of Drive Medical</small>	Hudson® Salter®
TUB NK 07	1112	2002
TUB NK 14	1118	2014
TUB NK 25	1119	2025
TUB NK 40	1124	2040
TUB NK 50	1120	2050



Universal Connector



Oxygen Tubing

## Oxygen Tubing with Universal Connector

TUB NK 07 U	7' Length, 50/cs
TUB NK 25 U	25' Length, 25/cs



Mask 002A

## Oxygen Mask

Transparent mask made with soft anatomical form for better fit. Wide adjustable elastic strap. Available in adult or pediatric sizes.

MASK 002A	Adult, 50/cs
MASK 002P	Pediatric, 50/cs

COMPARE TO:		
	MEDX <small>A Division of Drive Medical</small>	Hudson® Salter®
Mask 002A	1040	8115
Mask 002P	1035	1123



Mask 002P

## Oxygen Mask with 7' Tubing

Transparent oxygen mask is made with 7 feet of crush resistant safety tubing to prevent occlusion of oxygen supply. Comes on both adult and pediatric sizes.

MASK 003 A	Adult, 50/cs
MASK 003 P	Pediatric, 50/cs



Mask 003A

All of our tubing is made in translucent green to be in compliance with JACO/OSHA tripping hazard recommendations.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.



Ask for MEDX brand of Drive for a complete line of respiratory disposables.

respiratory

To order call toll free: 877.224.0946

drive  
WWW.DRIVEMEDICAL.COM