



# MABIS® NebPak™ UltraSonic Nebulizers

## Cool for School!

Introducing the first ultrasonic nebulizers designed specifically for kids, teens and young adults. NebPak Ultrasonic Nebulizers are available in bright colors with stylish camouflage bags, and are less than 8-oz., making these nebulizers ideal for school, sports or extracurricular activities.

- Small, quiet unit
- Handheld and lightweight
- Treatment time within 6 minutes for most medications
- Produces a fine mist for advanced medication therapy
- 3 options for nebulization rate
- Particle size range: 1 to 10 microns
- Automatic shut-off after 10 minutes
- Includes: mouthpiece, extension tube, connector, nasal piece, adult mask, child mask, medication cups, air filters, nebulizer stand, rechargeable battery pack, AC adapter, DC adapter and camouflage tote bag
- 5-year limited warranty

Item #	Description	Qty
40-365-000	 Green	1 ea
40-366-000	 Pink	1 ea



### NebPak Ultrasonic Nebulizer - Green

<b>Item Number</b>	40-365-000
<b>Unit UPC</b>	7-67056-40365-9
<b>Packaging Dimensions</b>	9-1/2" x 8" x 5"
<b>Packaging Number</b>	92-365-000
<b>Unit Weight</b>	Less than 8-oz.
<b>Unit Size</b>	2-3/8" x 4-5/8" x 6-1/2"

### NebPak UltraSonic Nebulizer - Pink

<b>Item Number</b>	40-366-000
<b>Unit UPC</b>	7-67056-40366-6
<b>Packaging Dimensions</b>	9-1/2" x 8" x 5"
<b>Packaging Number</b>	92-366-000
<b>Unit Weight</b>	Less than 8-oz.
<b>Unit Size</b>	2-3/8" x 4-5/8" x 6-1/2"